**BIOHIT**

Innovating for Health

Biohit Oyj

# **HELICOBACTER PYLORI IgG**

**INSTRUCTIONS FOR USE  
KÄYTTÖOHJE  
ANVÄNDNINGSSINSTRUKTIONER  
ARBEITSANLEITUNG  
MANUEL D'INSTRUCTION  
INSTRUCCIONES DE USO  
ISTRUZIONI PER L'USO  
ИНСТРУКЦИЯ ПОЛЬЗОВАТЕЛЯ**


Biohit Oyj  
Laippatie 1  
FIN-00880 Helsinki, Finland  
Tel: +358-9-773 861  
Fax: +358-9-773 86200  
E-mail: info@biohit.com  
www.biohit.com

601 040.02

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APPENDIX: QUALITY CONTROL CERTIFICATE

## 1. INTENDED USE

*Helicobacter pylori* IgG ELISA test is an enzyme-linked immunosorbent assay (ELISA) for the qualitative detection of human IgG class antibodies to *Helicobacter pylori* in EDTA plasma and serum. The test is intended to aid in the diagnosis of *Helicobacter pylori* infection in adult patients with clinical symptoms of gastritis. FOR *IN VITRO* DIAGNOSTIC USE.

## 2. SUMMARY AND EXPLANATION

*Helicobacter pylori* infection is the most important cause of chronic gastritis. Another mechanism for gastritis and severe atrophic gastritis is the autoimmune mechanism (1, 2). This kit is intended to aid in the diagnosis of *H. pylori* infection.

*Helicobacter pylori* is a spiral shaped, gram-negative bacterium that colonizes in the human stomach. The organism is found in the mucous layer of the stomach overlying the gastric epithelium and it does not appear to invade tissue. However, the mucosa underneath the area of the *H. pylori* colonization is invariably inflamed; this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life (1). This chronic inflammatory process can lead to atrophic gastritis, which has been linked with peptic ulceration and gastric cancer, two important diseases of the upper gastrointestinal tract (3-6). The presence of antibodies to CagA *H. pylori* strains have been linked with the development of atrophic gastritis in corpus (7). The epidemiological evidence indicates a link between gastric adenocarcinoma, mucosa associated lymphoid tissue (MALT) lymphoma (8,9) and *H. pylori* infection.

## 3. PRINCIPLE OF THE TEST

This *H. pylori* IgG antibody ELISA test is based on an enzyme immunoassay technique with partially purified *H. pylori* bacterial antigen adsorbed on a microplate and a detection antibody labeled with horseradish peroxidase (HRP).

The assay proceeds according to the following reactions:

1. Partially purified *H. pylori* bacterial antigen attached to the polystyrene surface of the wells binds *H. pylori* IgG antibodies present in the sample.

2. Wells are washed to remove residual sample.
3. HRP-conjugated monoclonal anti-human IgG binds to the *H. pylori* IgG antibodies.
4. The wells are washed and the TMB-substrate is added. The substrate is oxidized by the HRP enzyme, resulting in the formation of blue end product.
5. The enzyme reaction is terminated with the stop solution. *H. pylori* positive samples turn yellow with calculated values of >30 EIU.

#### **4. WARNINGS AND PRECAUTIONS**

**For *in vitro* diagnostic use.**

**CAUTION: Handle plasma and serum samples as potential biohazardous material.**

All samples should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the U. S. department of Health and Human Services (Bethesda, MD., USA) publication Biosafety in Microbiological and Biomedical Laboratories, 1999, 4th ed. (CDC/NIH) and No. (CDC) 88-8395 or any other local or national regulation on reports of laboratory safety procedures on different diseases.

This kit contains reagents manufactured from human blood components. The source materials provided in this kit have been tested for the presence of antibodies to hepatitis B and C, as well as antibodies to HIV, and found to be negative. However, as no test method can offer absolute assurance that these pathogens are absent, all recommended precautions for the handling of blood derivatives should be observed.

Always use protective gloves when handling patient samples. Use a safety pipetting device for all pipetting. Never pipette by mouth. Read all instructions prior to performing this assay. Dispose of all reagents provided in the kit by pouring into a sink and flushing with an excess of tap water.

#### **5. SPECIMEN COLLECTION AND HANDLING**

Blood sample is collected by venipuncture into e.g. a plastic EDTA or

serum tube. Blood tubes for plasma are mixed immediately by turning them upside down 5-6 times, and tubes for serum are allowed to clot (for minimum 30 minutes) at room temperature (20-25°C). Serum is separated after clotting, and plasma is separated immediately by centrifugation (e.g.) plastic tube, relative centrifugal force up to 2000 G, 10-15 minutes). Plasma/serum can be stored refrigerated (2-8 °C) for three days. For longer storage the samples should be stored frozen (preferably at -70°C, alternatively at -20°C). Mix the samples thoroughly after thawing. Avoid repeated freezing and thawing of the samples. Grossly hemolysed, lipemic or turbid specimen should be avoided.

## **6. KIT CONTENTS, REAGENT PREPARATION AND STABILITY OF THE MATERIALS PROVIDED**

The reagents supplied are sufficient for 96 wells and three separate runs. Reagents of different kit lots should not be mixed.

### **6.1. Microplate**

**Contents:** 12 x 8 strips in frame coated with partially purified *H. pylori* bacterial antigen.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date. Discard the strips after use.

### **6.2. Washing Buffer Concentrate (10x)**

**Contents:** 120 mL of 10 x phosphate buffered saline (PBS) concentrate containing Tween 20 and 0.1% ProClin 300 as preservative.

**Preparation:** Dilute 1 to 10 (e.g. 100 mL + 900 mL) with distilled water and mix well.

**Stability:** The diluted solution is stable for two weeks refrigerated (2-8°C).

### **6.3. Diluent Buffer**

**Contents:** 100 mL of phosphate buffer containing blocking protein, Tween 20, 0.1% ProClin 300 as preservative and red dye extract.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date.

### **6.4. Calibrator**

**Contents:** One vial containing 1.5 mL of human serum-based *H. pylori* IgG calibrator with 0.1% ProClin 300 as preservative. The EIU-value of the calibrator is printed on the label of the vial.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date.

### 6.5. Negative Control

**Contents:** One vial containing 1.5 mL of human serum-based *H. pylori* IgG negative control with 0.1% ProClin 300 as preservative. The EIU-value of the control serum is printed on the label of the vial.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.6. Positive Control

**Contents:** One vial containing 1.5 mL of human serum-based *H. pylori* positive control with 0.1% ProClin 300 as preservative. The EIU-value of the control serum is printed on the label of the vial.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date.

### 6.7. Conjugate Solution

**Contents:** 0.2 mL of HRP-conjugated monoclonal anti-human IgG in stabilizing buffer with 0.02% methylisothiazolone and 0.02% bromonitrodioxane, and 0.002% other active isothiazolones as preservatives.

**Preparation:** Dilute 1 to 100 (e.g. 40  $\mu$ L + 3960  $\mu$ L) with the diluent buffer.

**Stability:** Stable until expiry date. Discard the excess of the diluted conjugate solution after use.

### 6.8. Substrate Solution

**Contents:** 15 mL of tetramethylbenzidine (TMB) in aqueous solution.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date. Avoid exposure to direct light.

### 6.9. Stop Solution

**Contents:** 15 mL of 0.1 mol/L sulfuric acid.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date.

### 6.10. Incubation Covers

Four plastic sheets to cover the microplate during incubation.

### 6.11. Instructions for Use

## 7. MATERIALS REQUIRED BUT NOT PROVIDED

- ▶ Distilled or deionized water
- ▶ Micropipettes and disposable tips, to accurately deliver 5 – 1000  $\mu\text{L}$
- ▶ Pipettes to accurately deliver 1-10 mL
- ▶ 8-channel pipette delivering 100  $\mu\text{L}$
- ▶ Graduated cylinder, 1000 mL
- ▶ Vortex mixer for sample dilutions
- ▶ Test tubes for specimen dilutions
- ▶ Microplate washer
- ▶ Paper towels or absorbent paper
- ▶ Timer
- ▶ Incubator, 37°C
- ▶ Microplate reader, 450 nm
- ▶ e.g. plastic blood collection tube for serum or plasma
- ▶ Container for ice

## 8. STORAGE AND STABILITY

Store the *Helicobacter pylori* IgG kit refrigerated (2-8 °C). When stored at these temperatures, the kit is stable until the expiration date printed on the box label and the label of each individual kit component. Do not freeze or expose the kit to high temperatures, or store above 8°C when not in use. The substrate solution is sensitive to light. The microplate or individual strips should not be removed from the foil pouch until equilibrated to room temperature (20-25 °C). Unused strips must be returned to the foil pouch, sealed and stored at 2-8°C.

Do not use reagents after the expiration date printed on the label. Do not use reagents from kits with different lot numbers or substitute reagents from kits of other manufacturers. Use only distilled or deionized water. The components of the kit are provided at precise concentrations. Further dilution or other alterations of the reagents may cause incorrect results.

### Indication of Kit Deterioration

Liquid components should not be visibly cloudy or contain precipitated material. At 2-8°C, the washing buffer concentrate may, however, partially crystallize, but the crystals will dissolve by mixing at room temperature (20-25°C). The substrate solution should be colorless or pale blue. Any other color indicates deterioration of the substrate solution.

## Storage and stability of the samples

The results of the storage and freeze-thaw study showed that no differences in stability were observed in samples that were *H. pylori* negative, borderline, or clearly *H. pylori* positive. The stability of low positive specimens was not evaluated.

## 9. TEST PROCEDURE

### PRELIMINARY PREPARATIONS

Allow all reagents and the microplate to reach room temperature (20-25°C). Warm the incubator to 37°C. Dilute the washing buffer concentrate 1 to 10 (e.g. 100 mL + 900 mL) with distilled or deionized water. **Read the complete assay procedure before starting. It is recommended that the calibrator, controls and samples are applied to the plate in duplicate wells. It is necessary to use the calibrator and controls in each test run.**

Mix all the reagents well before use.

### STEP 1: SPECIMEN DILUTION

Dilute the mixed plasma or serum samples 1 to 200 (5 µL + 995 µL) with the diluent buffer, **mix well**.

### STEP 2: SAMPLE

Mix and pipette 100 µL of the diluent buffer (Blank), the calibrator (CAL), the negative control (NC), positive control (PC) and diluted samples (S1, S2, S3 etc.) into the wells as duplicates (See Figure 1). Cover the plate with the incubation cover. Incubate for 30 minutes at 37°C. Note: It is recommended that the samples are dispensed into the wells within 10 minutes to avoid assay drift within the plate.

	1	2	3	4	5
A	Blank	Blank			
B	CAL	CAL			
C	NC	NC			
D	PC	PC			
E	S1	S1			
F	S2	S2			
G	S3	S3			
H	etc.	etc.			

Figure 1. Pipetting Order.

### STEP 3: WASHING

Wash the strips with 3 x 350  $\mu\text{L}$  of the diluted washing buffer, and gently tap the inverted plate a few times on a clean paper towel to remove any residual washing buffer.

### STEP 4: CONJUGATE

Pipette 100  $\mu\text{L}$  of the mixed dilution (1 to 100) of conjugate solution into the wells with an 8-channel pipette. Cover the plate with the incubation cover. Incubate for 30 minutes at 37°C.

### STEP 5: WASHING

Wash the strips with 3 x 350  $\mu\text{L}$  of the diluted washing buffer, and gently tap the inverted plate a few times on a clean paper towel to remove any residual washing buffer.

### STEP 6: SUBSTRATE

Pipette 100  $\mu\text{L}$  of the substrate solution into the emptied wells with an 8-channel pipette. Start the incubation timer after pipetting the substrate solution into the first strip and continue the incubation for 30 minutes at room temperature (20-25°C). Avoid direct exposure to light during incubation.

### STEP 7: REACTION STOP

Pipette 100  $\mu\text{L}$  of the stop solution with an 8-channel pipette into the wells.

### STEP 8: MEASURING OF RESULTS

Measure the absorbance at 450 nm within 30 minutes.

## 10. RESULTS

### 10.1. Quality Control Values

1) The calibrator and controls included in the kit should be used in each test run. Quality control charts should be maintained to follow the performance of the controls or appropriate statistical methods should be used for analyzing control values, which should fall within the appropriate confidence intervals employed in each laboratory. The requirement is that positive control would be positive and negative control negative according to the values established in each laboratory. Expected control results must be obtained in order that the test results will be accepted.

2) Proper functioning of the test can be verified by observing the absorbance value of the calibrator. The absorbance reading of the calibrator should be > 1.000.

## 10.2. Calculation of the Results

Calculate the mean absorbance (A) value for the calibrator, controls and patient samples assayed in duplicate. Subtract the blank value from the calculated mean absorbance values.

Calculate the enzyme immunounits (EIU) as follows:

$$\frac{\bar{X}(A_{\text{Sample}}) - \bar{X}(A_{\text{Blank}})}{\bar{X}(A_{\text{Calibrator}}) - \bar{X}(A_{\text{Blank}})} \times 100 = \text{Sample EIU}$$

## 10.3. Interpretation of the Results

Negative < 30 EIU

Positive ≥ 30 EIU

A value less than 30 EIU indicates a negative result, or that no *H. pylori* antibodies have been detected, or that antibodies are below detectable level. A value of 30 EIU or over indicates that *H. pylori* antibodies are detected, and the result is thus positive. For quality purposes, negative and positive controls provided with the kit must be tested in each assay. Expected control results must be obtained in order to accept test results. The magnitude of the measured result above the cut-off value is not indicative of the total amount of antibody present. The cut-off values have been determined using the Biohit *Helicobacter pylori* IgG ELISA kit. From patients whose samples yield values near the cut-off, a second sample should be collected, if possible, within a reasonable time period. Each laboratory should establish its own range of expected values for the clinical situation where *H. pylori* IgG antibody values are used for diagnosis. Also the *H. pylori* IgG results determined for a given specimen or with assays from different manufacturers may vary due to differences in assay methods and reagent specificity. Results obtained by other assay methods from other manufacturers should not be used interchangeably.

## 10.4. Expected Values

Most persons exposed to *H. pylori* develop IgG antibodies to the organism (2,5,6). Age-specific rates for the presence of *H. pylori* antibodies are

similar for males and females. The prevalence of *H. pylori* infection is 30-40% in the U.S., Canada, and in Western Europe, about 20% in Australia, and 70-90% in Eastern Europe, Africa, South America and Asia. Many patients who have positive levels of antibody are asymptomatic, even though they are colonized with *H. pylori*. Therefore, antibody levels do not necessarily correlate with the severity of clinical symptoms.

## **11. LIMITATIONS OF THE PROCEDURE**

As with any diagnostic procedure the Biohit *Helicobacter pylori* IgG ELISA results must be interpreted together with the patient's clinical presentation and any other information available to the physician. A positive test for elevated antibodies may be achieved in patients with a recently treated *H. pylori* infection.

## **12. PERFORMANCE CHARACTERISTICS**

### **Reproducibility and Imprecision**

#### Inter-assay variation

A panel consisting of ten (10) samples (one negative, seven borderline, one positive and one high positive sample) was used for the determination of inter-assay variation. The samples were measured in duplicates in forty (40) different runs in 20 working days.

For serum samples, the mean coefficient of variation (%CV) ranged from 4.2 to 7.0% and for EDTA samples from 4.2 to 8.1%. The average %CV for borderline serum samples was 6.3% (range 5.4 – 7.0%) and for borderline EDTA samples 6.1% (range 5.0 – 8.1%).

#### Intra-assay variation

A panel consisting of four (4) serum samples (one negative, one borderline, one positive and one high positive) was used for the determination of intra-assay variation. The samples were measured in twenty two (22) single samples in one run and the average from three (3) different runs was calculated.

For serum samples, the mean coefficient of variation (%CV) ranged from 3.8 to 4.6% and for EDTA samples from 3.5 to 4.9%.

## Percent Agreements

The Biohit *Helicobacter pylori* IgG ELISA kit was evaluated in three studies with 523, 187 and 134 patient samples from age group of 20 to 89 years.

When compared to culture, the *Helicobacter pylori* IgG ELISA kit exhibited the following performance:

N=523		Culture	
		Positive	Negative
Biohit Diagnostics	Positive	166	9
	Negative	7	341

Clinical Sensitivity	96.5% (95%CI= 92.6 - 98.4%)
Clinical Specificity	96.0% (95%CI= 93.4 - 97.6%)

When compared to histology, the *Helicobacter pylori* IgG ELISA kit exhibited the following performance:

N=187		Histology	
		Positive	Negative
Biohit Diagnostics	Positive	52	8
	Negative	9	118

Percent Positive Agreement	85.2%
Percent Negative Agreement	93.7%

The Biohit *Helicobacter pylori* IgG ELISA test was compared to another commercial ELISA for *H. pylori* antibodies and exhibited the following performance:

N=134		Commercial ELISA	
		Positive	Negative
Biohit Diagnostics	Positive	32	1
	Negative	3	98

Percent Positive Agreement	91.4%
Percent Negative Agreement	99.0%

### **Analytical sensitivity**

The sensitivity of the *H. pylori* IgG ELISA was determined according to the NCCLS (CLIA) guideline Consensus Standards for Medical Testing (NCCLS EP17-A).

The limit of blank (LoB) was determined as the 95th percentile under which the measurement results of 60 blank samples (blank solution) fall. The LoB was found to be 1.1 EIU.

The limit of detection (LoD), which is the lowest actual amount of analyte that may be detected with 95% confidence, was assessed using four EDTA plasma and serum samples that contain antibody levels close to the observed LoB value. The LoD was found to be 1.8 EIU for both EDTA plasma and serum samples.

The limit of quantitation (LoQ), i.e. the lowest amount of analyte in a sample that can be reliably detected, was calculated assuming a coefficient of variation (CV%) of 20% at the level of LoD. Accordingly, the LoQ was calculated to be  $1.2 \times \text{LoD} = 2.2$  EIU for both EDTA plasma and serum samples.

## **13. DATE OF ISSUE**

### ***Helicobacter pylori* IgG ELISA kit insert.**

Version 01, 25.06.2009.

## **14. WARRANTY**

The Manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer. In case of interpretation disputes the English text applies.

All Biohit diagnostic kits have been manufactured according to our ISO 9001 / ISO 13485 quality management protocols and have passed all relevant Quality Assurance procedures related to these products.

## 15. SHORT OUTLINE OF THE PROCEDURE

Allow all the reagents to reach room temperature (20-25°C)  
Remember to mix all the reagents and samples well just before pipetting

\*

After mixing, pipette 100 µL of the blank solution, the calibrators, the control and diluted (1 to 200) patient samples into the wells

\*

Incubate for **30 min at 37°C**

\*

Wash the wells 3 times with 350 µL of the diluted washing buffer

\*

Pipette 100 µL of the diluted (1 to 100 dilution) and mixed conjugate solution into the wells

\*

Incubate for **30 min at 37°C**

\*

Wash the wells 3 times with 350 µl of the diluted washing buffer

\*

Pipette 100 µL of the mixed substrate solution into the wells

\*

Incubate for **30 min** at room temperature (**20-25°C**)

\*

Pipette 100 µL of the mixed stop solution into the wells

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Biohit Japan Co., Ltd.
Biohit Healthcare Limited
Biohit. Inc.