Biohit Oyj

Helicobacter pylori Quick Test

Quick Test Kit for the Detection of
H. pylori Urease Activity in
Biopsy Specimens

For *in vitro* diagnostic use
Store at 2-8°C Upon Receipt

**INSTRUCTIONS FOR USE
KÄYTTÖOHJE
ANVÄNDNINGSSINSTRUKTIONER
ARBEITSANLEITUNG
MANUEL D'INSTRUCTION
ISTRUZIONI PER L'USO
INSTRUCCIONES DE USO
ИНСТРУКЦИЯ ПОЛЬЗОВАТЕЛЯ**

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INSTRUCTIONS FOR USE***Helicobacter pylori* Quick Test****Cat. No. 602 015**

CONTENTS

1. INTENDED USE.....	6
2. CLINICAL BACKGROUND.....	6
3. PRINCIPLE OF THE QUICK TEST.....	6
4. WARNINGS AND PRECAUTIONS	7
5. SPECIMEN COLLECTION AND HANDLING	7
6. KIT CONTENTS AND REAGENT PREPARATION	7
6.1. Test plates	8
6.2. Instructions for Use	8
7. MATERIALS REQUIRED, BUT NOT PROVIDED.....	8
8. STORAGE AND STABILITY	8
9. TEST PROCEDURE	8
10. RESULTS	10
11. LIMITATIONS.....	10
12. PERFORMANCE CHARACTERISTICS.....	10
13. DATE OF ISSUE	11
14. WARRANTY	11
15. SHORT OUTLINE OF THE PROCEDURE	12

APPENDIX: QUALITY CONTROL CERTIFICATE

1. INTENDED USE

The *Helicobacter pylori* Quick Test kit is based on a biochemical reaction to qualitatively determine *H. pylori* urease activity directly from biopsy specimens at the laboratory or in Near Patient Testing (NPT) settings.

2. CLINICAL BACKGROUND

Helicobacter pylori infection is the most important cause of chronic gastritis. Another mechanism that may cause gastritis and severe atrophic gastritis is the autoimmune mechanism, which can also be triggered by an *H. pylori* infection (1, 2). This kit is intended to aid in the diagnosis of *H. pylori* infection.


Helicobacter pylori is a spiral shaped, gram-negative bacterium that colonizes the human stomach. The organism is found in the mucous layer of the stomach overlying the gastric epithelium, and it does not appear to invade tissue. However, the mucosa underneath the area of the *H. pylori* colonization is invariably inflamed; this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life (1). The chronic inflammatory process can lead to atrophic gastritis, which has been linked with peptic ulceration and gastric cancer, two of the most important diseases of the upper gastrointestinal tract (3-6).

The epidemiological evidence of a link between *H. pylori* infection and gastric adenocarcinoma or mucosa associated lymphoid tissue (MALT) lymphoma has resulted in classification of the organism as a group I carcinogen (7-9).

3. PRINCIPLE OF THE QUICK TEST

The quick test developed by Biohit for detecting *Helicobacter pylori* infection in stomach is based on the activity of the urease enzyme in a biopsy specimen. The biopsy specimen taken from the stomach is tested immediately. The development of the color in the test gel informs whether urease enzyme is present in the biopsy sample or not. In strong *H. pylori* colonization (high density) the reaction time is 2 minutes, in low density *H. pylori* colonization the reaction can take up to 30 minutes.

The assay proceeds according to the following reaction:

H. pylori contains a large amount of urease, which degrades urea to ammonia (NH₄⁺). The ammonia formed is detected by an indicator color present in the gel.

4. WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use.

CAUTION: Handle biopsy specimens as potential biohazardous material.

All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the U. S. department of Health and Human Services (Bethesda, MD., USA) publication Biosafety in Microbiological and Biomedical Laboratories, 1999, 4th ed. (CDC/NIH) and No. (CDC) 88-8395 on reports of laboratory safety procedures on different diseases or any other local or national regulation.

Always use protective gloves when handling patient samples. Read all instructions prior to performing the quick test. Do not use test plates beyond the expiry date. Discard the used test plates to biohazardous waste according to local and national regulations.

5. SPECIMEN COLLECTION AND HANDLING

Each test plate requires one biopsy specimen, which is recommended to be taken with 5 mm forceps either from the greater curvature of the middle antrum or the corpus. In order to assure high sensitivity, it is recommended that the *Helicobacter pylori* Quick Test is done by using biopsies from both antrum and corpus. (*H. pylori* colonization in stomach may sometimes be limited to either antrum or corpus). In performing the gastroscopy, the biopsies for the *Helicobacter pylori* Quick Test should be taken as early as possible to avoid possible errors caused by neutral and alkaline duodenogastric reflux. Before performing the test, blood should be removed from the biopsy specimen by placing it briefly on a sterile gauze pad. This is because a high amount of blood ($\geq 3.5 \times 10^{10}$ erythrocytes/l) can interfere with the test.

6. KIT CONTENTS AND REAGENT PREPARATION

Each Biohit *Helicobacter pylori* Quick Test kit contains reagents for 50 tests. Store the kit components refrigerated (2-8°C) until use.

6.1. Test Plates

Contents: Two bags each containing 25 test plates each (including one extra plate). Each plate has a well containing a yellow indicator gel for the biopsy specimen.

Preparation: Ready for use.

Stability: Stable until the expiration date. Discard the plates after use.

6.2. Instructions for Use

7. MATERIALS REQUIRED, BUT NOT PROVIDED

- Forceps
- Timer
- Sterile gauze pads
- Gloves

8. STORAGE AND STABILITY

Store the *H. pylori* Quick Test kit refrigerated (2-8°C). When stored at these temperatures the kit is stable until the expiration date printed on the box label and the label of test plate bags. Do not freeze or expose the kit to high temperatures or store at above 8°C when not in use. Do not use test plates after the expiration date printed on the label.

9. TEST PROCEDURE

Read the complete assay procedure before starting.

STEP 1

UREASE REACTION

1. Allow the test plate to reach room temperature for at least 30 minutes before use.
2. Open the label covering the well on the plate. Make sure that the color of the reagent gel is yellow. With the forceps place the biopsy specimen into the gel of the well. See figure 1.
3. Incubate for 1-2 minutes at room temperature (20-25°C).





Figure 1. Placing the Biopsy Specimen into the Well.

STEP 2

After 1-2 minutes of incubation, assess the change of color against the color chart given on the label. When the color of the reagent gel turns from yellow to intensive red, the test is positive. If there is no change in color, continue incubation until 30 minutes. If the color of the gel does not turn intensive red after 30 minutes total incubation, the test is negative. See the example below.

NEGATIVE

POSITIVE

Any color that appears (or disappears) after 30 minutes of incubation is not relevant and has no diagnostic value. The test plate must be discarded after 30 minutes

10. RESULTS

The test results are to be interpreted as follows:

Reagent gel: Yellow NEGATIVE (No *Helicobacter pylori* infection)

Reagent gel: Intensive red POSITIVE (*Helicobacter pylori* infection)

11. LIMITATIONS

False negative results

- can occur if *H. pylori* inhibiting antibiotics have been taken during the last 2-4 weeks prior to the examination.
- may occur if acid inhibiting drugs (PPI or H²-blockers) have been taken prior to the examination.

As with any diagnostic procedure, the Biohit *Helicobacter pylori* Quick Test results must be interpreted in the light of the patient's clinical presentation and any other information available to the physician.

12. PERFORMANCE CHARACTERISTICS

The golden standard *H. pylori* culture test was compared to the Biohit *Helicobacter pylori* Quick Test results as well as histological *H. pylori* detections in biopsy specimens. The culture test, considered as the golden standard, had 57 biopsies with *H. pylori* negative and 63 biopsies with *H. pylori* positive results. Results are shown in the table below including the specificities and sensitivities.

Culture	The Biohit <i>Helicobacter pylori</i> Quick Test		Histology	
	Negative	Positive	Negative	Positive
57 Negative	50	7	49	8
63 Positive	4	59	6	57
Specificity	87.7 %		86.0 %	
Sensitivity	93.7 %		90.5 %	

13. DATE OF ISSUE

Helicobacter pylori Quick Test kit insert.
Version 05, 25.06. 2009.

14. WARRANTY

The Manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer.

This Biohit Diagnostic kit has been manufactured according to our ISO 9001 / ISO 13485 quality management protocols and has passed all relevant Quality Assurance procedures related to this product.

15. SHORT OUTLINE OF THE PROCEDURE

Allow the test plate to reach room temperature (20-25°C) for at least 30 minutes before use

*

Place the biopsy specimen into the well

*

Incubate for 1-2 minutes at room temperature (20-25°C)

*

Check the color in the well

*

If there is no change in color continue incubation until 30 minutes at room temperature (20-25°C)

*

Interpret the result with the color chart on the test plate label

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1. Varis K, Sipponen P. Gastritis. In: Principles and Practice of Gastroenterology and Hepatology. Gitnick G (ed.). Appleton & Lange, Connecticut, 1994; 85-197.
2. Sipponen P. *Helicobacter pylori* gastritis-epidemiology. J Gastroenterol 1997; 32:273-277.
3. Sipponen P, Marshall BJ. Gastritis and gastric cancer. Western countries. Gastroenterol Clin North Am 2000; 29:579-592.
4. Wadström T. An update on *Helicobacter pylori* Current Opinion in Gastroenterology 1995; 11:69-75.
5. Northfield TC, Mendall M, Goggin PC. *Helicobacter pylori* Infection, Pathophysiology, Epidemiology and Management. Kluwer Academic Press; Dordrecht: 1994.
6. Sipponen P. Update on the Pathologic Approach to the Diagnosis of Gastritis, Gastric Atrophy, and *Helicobacter pylori* and its sequelae. J Clin Gastroenterol 2001; 32:96-202.
7. Parsonnet J, Friedman GD, Vandersteen DP, Chang Y, Vogelman JH, Orentreich N, Sibley RK. *Helicobacter pylori* infection and the risk of gastric carcinoma. N Engl J Med 1991; 325:1127-1131.
8. Parsonnet J, Hansen S, Rodriguez L, Gelb AB, Warnke RA, Jellum E, Orentreich N, Vogelman JH, Friedman GD. *Helicobacter pylori* infection and gastric lymphoma. N Engl J Med 1994; 330:1267-1271.
9. IARC Working Group on the evaluation of carcinogenic risks to humans. *Helicobacter pylori*. Schistomes, liver flukes, and *Helicobacter pylori*: views and expert opinions of an IARC Working group on the Evaluation of Carcinogenic Risks to Humans. Lyon: IARC, 1994: 177-240.

**ORDERING INFORMATION / TILAUS INFORMAATIO /
BESTÄLLNINGSPERAKTION / BESTELLINFORMATION/
INFORMATION CONCERNANT LA COMMANDE /
INFORMACIÓN DE PEDIDOS / INFORMAZIONI PER LE
ORDINAZIONI / ИНФОРМАЦИЯ ДЛЯ ЗАКАЗА**

Helicobacter pylori Quick Test kit

Cat. No. 602 015

Headquarters / Pääkonttori / Huvudkontor / Hauptsitz / Siège principal/
Oficina principal / Sede centrale / Головные офисы

BIOHIT OYJ

Laippatie 1

00880 Helsinki, Finland

Tel: +358-9-773 861

Fax: +358-9-773 86200

E-mail: info@biohit.com

www.biohit.com

CHINA

Biohit Biotech (Suzhou) Co. Ltd.

Tel: +86-21-6248 5589

Fax: +86-21-6248 7786

E-mail: info.china@biohit.com

FRANCE

Biohit SAS

Tel: +33-1-3088 4130

Fax: +33-1-3088 4102

E-mail: commercial.france@biohit.com

GERMANY

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Tel: +49-6003-828 20

Fax: +49-6003-828 222

E-mail: info@biohit.de

JAPAN

Biohit Japan Co., Ltd.

Tel: +81-3-5822 0021

Fax: +81-3-5822 0022

E-mail: sales@biohit.co.jp

RUSSIA

Biohit OOO Saint Petersburg

Tel: +7-812-327 5327

Fax: +7-812-327 5323

E-mail: main@biohit.ru

Biohit OOO, Moscow Office

Tel: +7-495-614 9550

Fax: +7-495-613 5577

Email: taras.pravdoljubenko@biohit.ru

U. K.

Biohit Healthcare Limited

Unit 1 Barton Hill Way

Torquay, Devon TQ2 8JG

Tel. +44-1803 315 900

Fax: +44-1803 315 530

E-mail: info@biohit.co.uk

U. S. A.

Biohit. Inc.

Tel: +1-732-922-4900

Fax: +1-732-922-0557

E-mail: diagnostics.usa@biohit.com

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Biohit Deutschland GmbH