Biohit Oyj

Lactose Intolerance Quick Test 10 and 25 Tests

Rapid Diagnosis of Lactose Intolerance,
Hypolactasia of the Small Intestine by Point
of Care (POC) Detection of Lactase Activity
from the Biopsy Specimen

For *in vitro* diagnostic use
Store at 2-8°C Upon Receipt

**INSTRUCTIONS FOR USE
KÄYTTÖOHJE
BRUKSANVISNING
ARBEITSANLEITUNG
MANUEL D'INSTRUCTION
ISTRUZIONI PER L'USO
INSTRUCCIONES DE USO
ИНСТРУКЦИЯ ПОЛЬЗОВАТЕЛЯ**

602 010 (25 tests)
602 012 (10 tests)

EXPLANATION OF THE SYMBOLS USED IN LABELS

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1. INTENDED USE.....	6
2. CLINICAL BACKGROUND.....	6
3. PRINCIPLE OF THE QUICK TEST.....	7
4. WARNINGS AND PRECAUTIONS.....	7
5. SPECIMEN COLLECTION AND HANDLING.....	8
6. KIT CONTENTS AND REAGENT PREPARATION.....	8
6.1. Test plates.....	9
6.2. Substrate Solution, SUBS.....	9
6.3. Chromogen Solution, CHRO.....	9
6.4. Signal Reaction Solution, SIGN.....	9
6.5. Instructions for Use.....	9
7. MATERIALS REQUIRED, BUT NOT PROVIDED.....	9
8. STORAGE AND STABILITY.....	9
9. TEST PROCEDURE.....	10
10. RESULTS.....	11
11. LIMITATIONS.....	12
12. PERFORMANCE CHARACTERISTICS.....	12
13. DATE OF ISSUE.....	12
14. WARRANTY.....	13
15. SHORT OUTLINE OF THE PROCEDURE.....	14

1. INTENDED USE

The Lactose Intolerance Quick Test kit is based on biochemical reactions which qualitatively determine lactase enzyme activity directly from biopsy specimens.

2. CLINICAL BACKGROUND

Lactose intolerance (hypolactasia) results in inability to digest lactose, the predominant sugar of milk. This inability results from a decrease of the lactase enzyme activity, which is produced in the small intestine. Lactase breaks down milk sugar (lactose) into monosaccharides (glucose and galactose), which can then be absorbed into the bloodstream.

Lactase activity begins to decrease after weaning. One-third of the adult population worldwide retain the ability to digest lactose. It is estimated that approximately 15-20% of Western and Northern Europeans, and 90% of Asians, Africans and Native Americans suffer from lactose intolerance. Often the symptoms remain undiagnosed, and, thus lactose intolerance may remain untreated for years.

Lactose intolerance may cause the following gastrointestinal symptoms after the intake of lactose: nausea, flatulence, abdominal bloating, cramps and pain, as well as diarrhea.

The patients seeking medical treatment for the symptoms caused by lactose intolerance are often referred to an upper intestinal endoscopy. However, lactose intolerance cannot be diagnosed on the basis of gastroscopy alone or the microscopical evaluation of biopsy specimen. The lack of the lactase enzyme must be determined biochemically from a homogenized biopsy specimen, or on the basis of a lactose load test, or with a breath hydrogen test. Most of these tests are uncomfortable for the patient, time-consuming, and expensive. Recently a DNA test has also been developed for screening purposes (1). The test is based on the finding that in the neighbourhood of the lactase gene there is an area where normal persons have cytosine/thymine or thymine/thymine, whereas lactose intolerant subjects have the mutation cytosine/cytosine. Therefore this test does not indicate mutation in lactase gene and is not functional either. Until recently, some of the above investigations have been the only means for diagnosing this disorder. Biohit's unique Point of Care (POC) -test in connection with gastroscopy enables quick diagnosis of the patients suffering from lactose intolerance.

3. PRINCIPLE OF THE QUICK TEST

The Lactose Intolerance Quick Test developed by Biohit for detecting lactose intolerance (hypolactasia of the small intestine) is based on the activity of the lactase enzyme in a biopsy specimen. The biopsy specimen taken from the mucosa of the upper part of the small intestine is examined immediately. The development of the color in the test liquid after 20 minutes informs whether or not the lactase enzyme is present in the biopsy specimen.

In the case of a normal biopsy specimen, the color develops as the lactase enzyme of the biopsy specimen breaks down the milk sugar added to the test buffer. However, if there is no or only a slight color reaction it can be concluded that the patient suffers from lactose intolerance.

4. WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use.

CAUTION: Handle biopsy specimens as potential biohazardous material.

All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the U. S. department of Health and Human Services (Bethesda, MD., USA) publication Biosafety in Microbiological and Biomedical Laboratories, 1999, 4th ed. (CDC/NIH) and No. (CDC) 88-8395 on reports of laboratory safety procedures on different diseases or any other local or national regulation. In the United States for research use only.

Always use protective gloves when handling patient samples. Read all instructions prior to performing the quick test. Do not use test plates beyond the expiry date. Discard the used test plates to biohazardous waste according to local and national regulations.

- Read all instructions before performing the test
- Do not interchange reagents between kit lots
- Do not use kit components beyond the expiration date.
- Gloves must be worn
- The chromogen solution, Bottle 2 contains 30% acetic acid, which is corrosive and can cause irritation. In a case of skin contact wash thoroughly with water. In a case of contact with eyes, wash for at least 15 minutes with water and seek for medical attention.

- If pipettes are used instead of the drop bottles, use a safety pipetting device for all pipetting. Never pipette by mouth.
- Reagents containing 0.002 % Thiomersal should be disposed of according to local and national regulations.

5. SPECIMEN COLLECTION AND HANDLING

Biopsy specimens are recommended to be taken with e.g. 5 mm forceps from the mucosa of the upper part of the small intestine at any site of the postbulbar duodenal wall during gastroscopy. The bulbus should be avoided, because of the secondary changes in the lactase activity in the bulbar mucosa due to acid related diseases in some patients.

Total lactase content in the biopsy depends also to some extent on the size of the biopsy specimen used. Estimation of the recommended biopsy size can be done visually by comparing the size of each biopsy to the picture below.

Recommended size of biopsies:

Diameter (\varnothing) between

The performance of the Biohit Lactose Intolerance Quick Test was evaluated using biopsies in this size range. If the biopsy is smaller than recommended, there is a risk of a false negative (hypolactasia) result. If the biopsy specimen is too large, there is a risk of a false positive (normal) result.

Before performing the test, blood can be removed from the biopsy specimen by placing it briefly on the sterile gauze pad. It is recommended to perform the test immediately.

6. KIT CONTENTS AND REAGENT PREPARATION

Each Biohit Lactose Intolerance Quick Test kit contains reagents for 10 (Cat. no. 601 012) or 25 (Cat. no. 602 010) tests. Store the kit components refrigerated (2-8°C) until use.

6.1. Test Plates

Contents: 10 (Cat. no. 601 012) or 25 (Cat. no. 602 010) test plates in a bag, each containing a well for the biopsy specimen.

Preparation: Ready for use.

6.2. Substrate Solution, Bottle 1, SUBS

Contents: 1 vial containing 3.5 ml of acetate buffer and lactose with 0.002% Thiomersal as preservative.

Preparation: Ready for use.

Stability: Stable until expiration date.

6.3. Chromogen Solution, Bottle 2, CHRO

Contents: 1 vial containing 0.75 ml of chromogen solution in acetic acid with 0.002 % Thiomersal as preservative.

Preparation: Ready for use.

Stability: Stable until expiration date. **WARNING:** Acetic acid is corrosive and can cause irritation!

6.4. Signal Reaction Solution, Bottle 3, SIGN

Contents: 1 vial containing 3.5 ml of enzyme solutions with 0.002 % Thiomersal as preservative.

Preparation: Ready for use.

Stability: Stable until expiration date.

6.5. Instructions for Use

7. MATERIALS REQUIRED, BUT NOT PROVIDED

- Forceps
- Timer
- Sterile gauze pads
- Gloves

8. STORAGE AND STABILITY

Store the Lactose Intolerance Quick Test kit refrigerated (2-8 °C). When stored at these temperatures the kit is stable until the expiration date printed on the box label and the label of each individual kit component. Do not freeze or expose the kit to high temperatures or store at above 8°C when not in use.

Do not use reagents after the expiration date printed on the label. Do not use reagents from kits with different lot numbers or substitute reagents from kits of other manufacturers.

9. TEST PROCEDURE

Total reaction time is 20 minutes and the test is done in two steps: First a 15 minute lactase reaction followed immediately by a 5 minute signal reaction. Gloves must be worn. Allow all the reagents to reach room temperature (20-25°C) for at least 15 minutes before use. Mix all the reagents before using by turning the bottles upside down for a few times. Then tap the bottles on the table to ensure the liquid returns to the bottle. When adding the drops into the wells, hold the bottles in vertical upside down position (See figures 1-2).

STEP I

LACTASE REACTION

1. Open the label covering the well on the plate. Place the biopsy specimen with the forceps into the well. See figure 1.
2. Add 2 drops (80 µl) of the substrate solution (Bottle 1) into the well. Close the label.
3. Mix the plate by shaking 5-6 times vigorously sideways on the table.
4. Incubate for 15 minutes at room temperature (20-25 °C).


Figure 1. Lactase Reaction.

STEP II

SIGNAL REACTION

1. Open the label again. Add 1 drop (10 μ l) of the chromogen solution (Bottle 2) into the well. See figure 2.
2. Immediately add 2 drops (80 μ l) of the signal reaction solution (Bottle 3). Close the label.
3. Mix the plate by shaking for 5-6 times vigorously sideways on the table.
4. Incubate totally for 5 minutes at room temperature (20-25 $^{\circ}$ C) and interpret the lactase content with the color chart on the plate label immediately.


Figure 2. Signal Reaction.

10. RESULTS

The lactase content of the specimen can be interpreted with the color chart given on the label. See the example below.

-  Severe hypolactasia
-  Mild hypolactasia
-  Normolactasia

In case of normolactasia (lactase activity > 10 U/g protein and lactase/saccharase ratio > 0.25, U= μ mol/min, 37 $^{\circ}$ C) hypolactasia can be excluded. For the more detailed interpretation of the results please refer to the interpretation chart included in the kit.

11. LIMITATIONS OF THE PROCEDURE

Test results may be false if a biopsy of an incorrect size is used. The diameter of biopsies to be tested should be 1.5-2.0 mm.

As with any diagnostic procedure the Biohit Lactose Intolerance Quick Test results must be interpreted in the light of the patient's clinical presentation and any other information available to the physician.

12. PERFORMANCE CHARACTERISTICS

Correlation With Other Methods

The Biohit Lactose Intolerance Quick Test results of biopsy specimens were compared with the results of biopsy specimens in the biochemical determination of lactase activity in which >10 U/g indicates normolactasia and <10 U/g mild or severe hypolactasia (2).

Sensitivity: 95 %

Specificity: 100 %

Positive predictive value (PPV): 100 %

Negative predictive value (NPV): 98 %

27.5% of the patients had severe hypolactasia. The validity of the lactase activity was confirmed by determining the lactase/saccharase ratio. The number of symptomatic patients tested was 80.

When Lactose Intolerance Quick Test was compared to the combination of H₂ and CH₄ breath test results, the following characteristics were shown (3): sensitivity 96%, specificity 96%, positive predictive value (PPV) 96%, negative predictive value 96%, N=50.

13. DATE OF ISSUE

Lactose Intolerance Quick Test kit insert.

Version 01, 18.05.2010.

14. WARRANTY

The Manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer.

This Biohit Diagnostic kit has been manufactured according to our ISO 9001 / ISO 13485 quality management protocols and has passed all relevant Quality Assurance procedures related to this product.

15. SHORT OUTLINE OF THE PROCEDURE

Allow all the reagents to reach room temperature (20-25°C) for at least 15 minutes before use. Mix all the reagents before using by turning the bottles upside down for a few times.

*

Tap the bottles on the table to ensure the liquid returns to the bottle.

When adding the drops into the wells, hold the bottles in vertical upside down position.

*

Place the biopsy specimen into the well.

*

Add 2 drops (80 μ l) of the Substrate Solution (Bottle 1, SUBS) into the well. Close the label.

*

Mix the plate by shaking for 5-6 times vigorously sideways on the table.

*

Incubate for 15 minutes at room temperature (20-25°C).

*

Add 1 drop (10 μ l) of the Chromogen Solution (Bottle 2, CHRO) into the well and immediately 2 drops (80 μ l) of the Signal Reaction Solution (Bottle 3, SIGN). Close the label. Mix the plate by shaking 5-6 times vigorously sideways on the table.

*

Incubate for 5 minutes at room temperature (20-25°C).

*

Interpret the lactase content with the color chart on the label.

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Biohit Healthcare Limited
Biohit SAS
Biohit Deutschland GmbH
Biohit Japan Co., Ltd.