

BIOHIT

Innovating for Health

Biohit Oyj

**Lactase
CONTROL+**

Control Solution for
Lactose Intolerance
Quick Test



INSTRUCTIONS FOR USE

Biohit Oyj
Laippatie 1
FIN-00880 Helsinki, Finland
Tel: +358-9-773 861
Fax: +358-9-773 86200
E-mail: info@biohit.com
www.biohit.com

REF 602 018

IVD



Lactase Control+

Cat. No. 602 018

CONTENTS

1. INTENDED USE.....	3
2. WARNINGS AND PRECAUTIONS	3
3. PRODUCT CONTENTS AND STABILITY	3
4. USE OF THE LACTASE CONTROL+	3
5. DATE OF ISSUE	4
6. WARRANTY	5
7. EXPLANATION OF THE SYMBOLS USED IN LABELS.....	5
8. ORDERING INFORMATION	6

1. INTENDED USE

The use of the Biohit Lactase Control+ solution is intended as a quality control reagent for the Biohit Lactose Intolerance Quick Test (Cat. No. 602 010 and 602 012).

2. WARNINGS AND PRECAUTIONS

- Read all instructions before using Lactase Control+.
- Do not use Lactase Control+ beyond the expiration date.
- Gloves must be worn.
- If pipettes are used instead of the drop bottles, use a safety pipetting device for all pipetting. Never pipette by mouth.
- The solution contains 0.002 % Thiomersal, and it should be disposed of according to local and national regulations.

3. PRODUCT CONTENTS AND STABILITY

Contents: 1 drop bottle containing 1.0 ml of lactase solution with 0.002% Thiomersal as preservative.

Preparation: Ready for use.

Stability: When stored as directed (at 2-8 °C), Biohit Lactase Control+ will remain stable until the labeled expiration date.

4. USE OF THE LACTASE CONTROL+

The test procedure with the Lactase Control+ is similar to the Lactose Intolerance Quick Test reaction, except the Lactase Control+ is used instead of a biopsy sample. Please read the instructions carefully before starting the test procedure. Mix all the reagents before use by turning the bottles upside down a few times. Tap the bottles onto a table to ensure the liquid returns to the bottle. When adding the drops into the wells, hold the bottles in a vertical upside down position.

STEP I: LACTASE REACTION

1. Open the label covering the well on the plate. Add one drop (10 µl) of the Lactase Control+ solution into the well.
2. Add 2 drops (80 µl) of the substrate solution (Bottle 1) into the well. Close the label.
3. Turn the plate upside down twice by holding your thumb tight-

ly on top of the well.

4. Incubate for 15 minutes at room temperature (20-25 °C).

STEP II: SIGNAL REACTION

1. Open the label again. Add 1 drop (10 µl) of the chromogen solution (Bottle 2) into the well.
2. Immediately add 2 drops (80 µl) of the signal reaction solution (Bottle 3). Close the label.
3. Turn the plate upside twice by holding your thumb tightly on top of the well.
4. Incubate for 5 minutes at room temperature (20-25 °C).

STEP III: RESULTS

After 5 minutes of the signal reaction, the result can be interpreted with the color chart also given in the Lactose Intolerance Quick Test.

Valid: The color of the solution is clearly blue (Mild or Normolactasia in the color chart).

Invalid: There is no color or the color of the solution is only slightly blue (Hypolactasia in the color chart).



If the Lactase Control+ does not perform as expected, repeat the test or contact Biohit for technical assistance before testing patient samples. DO NOT use reagents beyond their labeled expiration dates.

5. DATE OF ISSUE

Lactase Control+ kit insert
Version 01, October 28, 2008.





6. WARRANTY

The Manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer. In case of interpretation disputes the English text applies.

The product has been manufactured according to our ISO 9001 / ISO 13485 quality management protocols and have passed all relevant Quality Assurance procedures related to these products.

7. EXPLANATION OF THE SYMBOLS USED IN LABELS

	Catalogue number
	Batch code
	Use by
 2-8 °C	Temperature limitation. Store at 2-8 °C

8. ORDERING INFORMATION

Headquarters
BIOHIT OYJ
Laippatie 1
00880 Helsinki, Finland
Tel: +358-9-773 861
Fax: +358-9-773 86200
Email: info@biohit.com
www.biohit.com

CHINA
Biohit Biotech (Suzhou) Co. Ltd.
Tel: +86-21-6248 5589
Fax: +86-21-6248 7786
E-mail: info.china@biohit.com

FRANCE
Biohit SAS
Tel: +33-1-3088 4130
Fax: +33-1-3088 4102
Email: commercial.france@biohit.com

GERMANY
Biohit Deutschland GmbH
Tel. +49-6003-828 20
Fax: +49-6003-828 222
E-mail: info@biohit.de

JAPAN
Biohit Japan Co., Ltd.
Tel: +81-3-5822 0021
Fax: +81-3-5822 0022
E-mail: sales@biohit.co.jp

RUSSIA
Biohit OOO Saint Petersburg
Tel: +7-812-327 5327
Fax: +7-812-327 5323
E-mail: main@biohit.ru

Biohit OOO, Moscow
Tel: +7-495-214 9550
Fax: +7-495-213 5577
Email: taras.pravdoljubenko@biohit.ru

U. K.
Biohit Ltd.
Tel. +44-1803 315 900
Fax: +44-1803 315 530
E-mail: sales@biohit.co.uk

U. S. A.
Biohit. Inc.
Tel: +1-732-922-4900
Toll Free: 800-922 0784
Fax: +1-732-922-0557
E-mail: diagnostics.usa@biohit.com