# BIOHIT HealthCare

Innovating for Health

## HELICOBACTER PYLORI quick test

for the detection of H.pylori in biopsy specimens

REF 602015 (50 tests)

## ENGLISH

### 1. INTENDED USE

The Helicobacter pylori quick test kit is based on a biochemical reaction to qualitatively determine H. pylori urease activity directly from biopsy specimens in Near Patient Testing (NPT) settings (eg. Endoscopy) and laboratories.

### 2. CLINICAL BACKGROUND

Helicobacter pylori infection is the most important cause of chronic gastritis. Another mechanism that may cause gastritis and severe atrophic gastritis is the autoimmune mechanism, which can also be triggered by an H. pylori infection [1, 2]. This kit is intended for the diagnosis of H. pylori infection. Helicobacter pylori is a spiral shaped, gram-negative bacterium that colonizes the human stomach. The organism is found in the mucous layer of the stomach overlying the gastric epithelium. The mucosa colonized by H. pylori is invariably inflamed; this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life [1]. The chronic inflammatory process can lead to atrophic gastritis, which has been linked with peptic ulceration and gastric cancer. Two of the most important diseases of the upper gastrointestinal tract [3-6]. The epidemiological evidence of a link between H. pylori infection and gastric adenocarcinoma or mucosa associated lymphoid tissue (MALT) lymphoma has resulted in classification of the organism as a group I carcinogen [7-9].

### 3. PRINCIPLE OF THE QUICK TEST

The quick test developed by Biohit for detection of Helicobacter pylori infection in stomach is based on the activity of the urease enzyme in a biopsy specimen. The biopsy specimen taken from the stomach is tested immediately. The development of the color in the test gel informs whether urease enzyme is present in the biopsy sample or not. In strong H. pylori colonization (high density) the reaction time is 2 minutes, in low density H. pylori colonization the reaction can take up to 30 minutes. The assay proceeds according to the following reaction: (NH<sub>4</sub>)<sub>2</sub>CO<sub>3</sub> + 2H<sub>2</sub>O + H<sup>+</sup> →→ 2NH<sub>4</sub><sup>+</sup> + HCO<sub>3</sub><sup>-</sup>. H. pylori produces urease, which degrades urea to ammonia (NH<sub>3</sub>). The ammonia formed is detected by an indicator color present in the gel.

### 4. WARNINGS AND PRECAUTIONS

For in vitro diagnostic use. Handle all specimens as potential biohazardous material. All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the U. S. department of Health and Human Services (Bethesda, MD, USA) publication Biosafety in Microbiological and Biomedical Laboratories, 1999, 4th ed. [CDC/NIH] and No. (CDC) 88-8395 on reports of laboratory safety procedures on different diseases or any other local or national regulation. Always use protective gloves when handling patient samples. Read all instructions prior to performing the quick test. Do not use test plates beyond the expiry date. Discard the used test plates to biohazardous waste according to local and national regulations.

### 5. SPECIMEN COLLECTION AND HANDLING

Each test plate requires one biopsy specimen, which is recommended to be taken with 5 mm forceps either from the greater curvature of the middle antrum or the corpus. In order to assure high sensitivity, it is recommended that the Helicobacter pylori quick test is done by using biopsies from both antrum and corpus. [H. pylori colonization in stomach may sometimes be limited to either antrum or corpus]. Several biopsies of one patient can be taken simultaneously on one plate, if they fit well into the gel. The biopsies for the Helicobacter pylori quick test (HPQT) should be taken as early as possible during gastroscopy, to avoid possible errors caused by neutral and alkaline duodenogastric reflux. Do not contaminate the biopsy forceps or biopsy specimen with formalin or other tissue fixatives. It is recommended that the biopsies for HPQT test are collected first, before any biopsies are collected for histology. If the forceps become contaminated with formalin or other tissue fixatives rinse thoroughly in water before collecting the biopsies for HPQT test. Before performing the test, blood should be removed from the biopsy specimen by placing it briefly on a sterile gauze pad. This is because a high amount of blood [≥ 3.5 x 10<sup>10</sup> erythrocytes/l] can interfere with the test.

### 6. KIT CONTENTS AND REAGENT PREPARATION

Each Biohit Helicobacter pylori quick test kit contains reagents for 50 tests. Store the kit components refrigerated (2-8°C) until use.

### 6.1. Test Plates

Contents: Two bags each containing 25 test plates each. Each plate has a well containing a yellow indicator gel for the biopsy specimen.

Preparation: Ready for use.

Stability: Stable until the expiration date.

### 6.2. Instructions for Use

### 7. MATERIALS REQUIRED, BUT NOT PROVIDED

- Forceps
- Timer
- Sterile gauze pads
- Gloves

### 8. STORAGE AND STABILITY

Store the H. pylori quick test kit refrigerated (2-8°C). When stored at this temperature the kit is stable until the expiration date printed on the box label and the label of test plate bags. Do not freeze or expose the kitto high temperatures or store at above 8°C when not in use. Do not use test plates after the expiration date printed on the label.

### 9. TEST PROCEDURE

Read the complete assay procedure before starting.

**STEP 1:** Urease reaction

- Allow the test plate to reach room temperature for at least 30 minutes before use.
- Open the label covering the well on the plate. Make sure that the color of the reagent gel is yellow. Place the biopsy into the gel of the well so that it is totally surrounded by the gel. See figure 1.
- Incubate for 1-2 minutes at room temperature [20-25°C].

### STEP 2

After 1-2 minutes of incubation, assess the change of color against the color chart given on the label. When the color of the reagent gel turns from yellow to intensive red, the test is positive. Any RED colour in the gel (except blood) indicates a POSITIVE result - lower density colonisations may not change the entire gel to RED. If there is no change in color after 2 minute, continue incubation for a maximum of 30 minutes. If the color of the gel does not turn intensive red after 30 minutes total incubation, the test is negative. See figure 2.

Any color that appears (or disappears) after 30 minutes of incubation is not relevant and has no diagnostic value. The test plate must be discarded after 30 minutes.

### 10. RESULTS

The test results are to be interpreted as follows:

Reagent gel: Yellow NEGATIVE (No Helicobacter pylori infection)

Reagent gel: Intensive red POSITIVE (Helicobacter pylori infection)

### 11. LIMITATIONS

False negative results

- can occur if H. pylori inhibiting antibiotics have been taken 2-4 weeks prior to the examination.
- may occur if acid inhibiting drugs (PPI or H2-blockers) have been taken prior to the examination.

As with any diagnostic procedure, the Biohit Helicobacter pylori Quick Test results must be interpreted in the light of the patient's clinical presentation and any other information available to the physician.

### 12. WARRANTY

Biohit shall remedy all defects discovered in any product [the "Defective Product"] that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in Biohit's specifications for the Products. ANY WARRANTIES WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS. CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL. The period of this warranty is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by Biohit. This Biohit Diagnostic kit has been manufactured according to ISO 9001/ ISO 13485 quality management protocols. In case of interpretation disputes the English text applies.

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