BIOHIT HealthCare

Innovating for Health

ENGLISH

BIOHIT *Helicobacter pylori* UFT300 quick test for the detection of *H. pylori* in biopsy specimens
REF 602 005 [5 tests], 2,5 ml of reagent solution
REF 602 019 [50 tests], 6,5 ml of reagent solution
REF 602 021 [100 tests], 13 ml of reagent solution
REF 602 005PLA [5 tests with plate], 2,5 ml of reagent solution
REF 602 019PLA [50 tests with plate], 6,5 ml of reagent solution

1. INTENDED USE AND BACKGROUND

The Biohit *Helicobacter pylori* UFT300 quick test aids the rapid diagnosis of *H. pylori* from human gastric mucosa. The qualitative *in vitro* test is based on the detection of bacterial urease activity in gastric biopsy specimens.

Helicobacter pylori infection is the most important cause of chronic gastritis. A gastric mucosa colonized by *H. pylori* is invariably inflamed, this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life [1]. The chronic inflammatory process can lead to atrophic gastritis, which has been linked with peptic ulcer and gastric cancer, two of the most important diseases of the upper gastrointestinal tract [2-5].

2. PRINCIPLE OF THE QUICK TEST

H. pylori is characterized by its capacity to produce urease, an enzyme able to break down urea into ammonium and bicarbonate, with the reaction: [NH₂]CO + 2H₂O + H⁺ ->> 2NH₄⁺ + HCO₃⁻. To selectively identify *H. pylori* with high sensitivity, the test exploits the pH variation caused by the production of ammonium, revealed by the color change of the pH indicator from yellow to magenta.

3. WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
2. The reagent solution bottle must be kept closed until the time of use.
3. All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious.
4. Contact of the test solution with skin or eyes should be prevented. If the solution contacts the skin or eye, flush with large amount of water.
5. In order to obtain accurate results, package insert instructions should be followed.
6. Do not use beyond the expiration date printed on the outside of the box.
7. The stability of the opened reagent bottle is 3 months.
8. Disinfection of the bottle is recommended to be done by spraying the reagent bottle cap closed and drying any residual disinfectant with a clean cloth from the surface of the bottle. Disinfectants (e.g. Deconex, Oxivir, ethanol) may cause false negative or false positive results if in contact with the reagent.
9. Dispose materials according to local waste management legislaion.

4. SPECIMEN COLLECTION AND HANDLING

Patients should not have taken antibiotics or bismuth salts within three weeks, or H2-blockers and PPI drugs within two weeks, preceding endoscopy. Incomplete eradication of *H. pylori* can give false negative results. It is recommended that biopsy specimens are tested immediately after collection. Do not let the specimen dry before testing. If there is visible blood or bile on the specimens, dab the biopsy briefly on sterile gauze. Even a small amount of blood can mask the color of the reaction.

Do not contaminate the biopsy forceps or biopsy specimen with lidocaine hydrochloride or formalin, or other lubricants or tissue fixatives. It is recommended that the biopsies for UFT300 test are collected first, before any biopsies are collected for histology. If the forceps become contaminated with formalin or other tissue fixatives rinse thoroughly in water before collecting the biopsies for UFT300 test. Contamination with lidocaine hydrochloride (e.g. Xylocaine®) can inhibit the reaction and cause false negative results.

5. TEST PROCEDURE

- Write the ID of the patient on the test tube’s label (provided on separate sheets)/ on plate label
- Open the tube/label covering the well of the plate and add one or more biopsies from the same patient into the tube/well.
- Add one drop of the reagent solution into the test tube/well – just the minimum amount required to cover the biopsy (Figure 1). Keep the bottle upright while dropping (Figure 2). Only if the biopsy is not fully covered by the solution, add one more drop; avoid adding too much of the reagent solution, do not to dilute the reaction. Put the cap back on the tube/put the label back to cover the well of the plate.
- Shake the tube/plate by swirling it horizontally on a flat surface for 5 seconds. Make sure the biopsy and the solution are at the bottom of the test tube/well of the plate by giving the tube/plate a light tap.

6. RESULTS

The color reaction is easy to interpret by removing the cap/label of the plate and looking into the tube from the top. The presence of *H. pylori* in the biopsy makes the color of the reagent solution turn from yellow to magenta. The presence of blood or alkaline substances can induce a slight variation of the intensity of yellow color, without causing a turn to magenta. The color should be read at 5 min from the biopsy insertion for both positive and negative results. Do not read the color after 5 min. Discard the tube/plate after 5 minutes reading. (See also invalid result below).

Negative result: Color of test solution remains yellow. Figure 3.
Positive result: Color of test solution turns from yellow to magenta. Figure 4.
Invalid result: If the color of the test solution turns from yellow to magenta immediately [within seconds] from the insertion of the biopsy, the result can be influenced by the presence of blood or alkaline substances (eg. bile) in the biopsy. It is suggested to add 1 or more drops of the reagent solution, until the starting yellow color is obtained, and then to shake well the tube, taking care that the biopsies remain fully immersed.

The time needed for the color change is inversely proportional to the *H. pylori* concentration in the biopsy. Analytical sensitivity of UFT300 was evaluated with different concentrations of *H. pylori* suspension: 10⁴ - 10⁸ bacterial cells of *H. pylori* turn the color to magenta at 5 minutes.

7. LIMITATIONS

- UFT300 test is a qualitative test for the detection of *H. pylori* from gastric biopsies.
- The result of the test - as for all diagnostic tests - must be considered together with other clinical data available to the physician.
- False negative results can occur when the concentration of *H. pylori* in the biopsy is below the sensitivity limit of the test, or when the biopsy does not contain *H. pylori*, because of the uneven distribution of bacteria in the gastric mucosa. In case of disagreement with other diagnostic criteria, the follow-up with histologic or cultural methods is recommended.
- Therapy with antibiotics, H2-blockers or proton pump inhibitors can adversely influence the test results.

8. WARRANTY

Biohit shall remedy all defects discovered in any Product (lhe "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in Biohit’s specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL. The period of this warranty is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by Biohit. This Biohit Diagnostic kit has been manufactured according to ISO 9001 / ISO 13485 quality management protocols. In case of interpretation disputes the English text applies.

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