Hyvä Kollega

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At an early stage	The stochastic GastroSoft report states:	¹³ C - urea breath test or Stool antigen testreport:
The diagnosis for		
Functional vs. organic dyspepsia.	YES	NO
When GastroPanel indicates the gastric mucosa is healthy, the dyspepsia complaints are often caused by functional dyspepsia or another disease not involving the gastric mucosa		
H. pylori infection (gastritis)	YES	NOT RELIABLE (1)
Atrophic gastritis (damaged and severely dysfunctional gastric mucosa) and the probabilities of different conditions affecting the mucosa of the gastric corpus or antrum or both (normal, gastritis or atrophic gastritis)	YES	NO
The risks (related to atrophic gastritis) of		
Gastric cancer	YES	YES/NO (2)
Vitamin B12 deficiency	YES	NO
Peptic ulcer disease	YES	YES/NO (3)
The risks of the complications of		
Gastroesophageal reflux disease:		
Esophagitis and Barrett's esophagus	YES (4)	NO
If necessary, a recommendation for		
Gastroscopy and biopsy examination	YES	NO
Treatment of <i>H. pylori</i> infection	YES	YES/NO (5)
Determination of vitamin B12 and homocysteine	YES	NO
Follow-up examination to monitor		
the incidence of atrophic gastritis	YES	NO
the healing of the <i>H. pylori</i> infection	YES	YES
the healing of atrophic gastritis	YES	NO

(1) The ¹³C- urea breath - and stool antigen tests give **false negative results** if the patient has **a) atrophic gastritis** (a risk of gastric cancer and peptic ulcer disease and vitamin B12 deficiency and related diseases, such as dementia, depression and polyneuropathies as well as atherosclerosis, strokes and heart attacks) or **b) MALT lymphoma** or **c) bleeding peptic ulcer disease** or **d) if the patient is currently receiving antibiotics or PPIs** (proton pump inhibitors).

(2) The risk of gastric cancer is very low without atrophic gastritis in corpus, antrum or both. But in some cases, a *H. pylori* infection without histologically observable atrophic gastritis may be associated with gastric cancer and peptic ulcer disease.

(3) No peptic ulcer disease with corpus atrophy (no acid, no ulcer). The risk of peptic ulcer disease is very low without antrum atrophy.

(4) High pepsinogen I (over 120 ug / l) and high pepsinogen I and II ratio I (over 10) and low gastrin-17 (below 2 umol / l) indicate high acid (HCl) output and risks for the complications of gastroesophageal reflux disease.

(5) When the incidence of *H. pylori*-related atrophic gastritis is monitored, the patient can be offered targeted, safe treatment at the right time. The need for medication and the costs and adverse effects of medication can thus be reduced. If the patient has been diagnosed with peptic ulcer disease (gastric or duodenal ulcer), the *H. pylori* infection has to be treated (6). It should also be treated if the patient has atrophic gastritis. The patient and the doctor may also agree on eradication treatment for other reasons for example when the patient's close relatives have been diagnosed with gastric cancer.

(6) Press Release: The 2005 Nobel Prize in Physiology or Medicine, 3 October 2005 jointly to Barry Marshall and J. Robin Warren for their discovery of "the bacterium *Helicobacter pylori* and its role in gastritis and peptic ulcer disease": - "An indiscriminate use of antibiotics to eradicate *Helicobacter pylori* also from healthy carriers would lead to severe problems with bacterial resistance against these important drugs. Therefore, treatment against *Helicobacter pylori* should be used restrictively in patients without documented gastric or duodenal ulcer disease." <http://nobelprize.org/medicine/laureates/2005/press.html>

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The European Helicobacter Pylori Study Group

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Viitteet:

1. Färkkilä M, Miten dyspepsia tulisi hoitaa, Duodecim 2004; 120: 2537– 42
2. Färkkilä M, Kaksiportainen menetelmä seuloa dyspepsia sairastavat, Wyeth News to Doctors 2005
3. In the two-stage screening method, the 13C urea breath test is recommended as the primary test for detecting Helicobacter pylori infections. This test is nevertheless only reliable when Helicobacter pylori infection is fresh and the number of bacteria is high. Its sensitivity decreases significantly if the number of bacteria is small. For this reason, the 13C urea breath test gives false negative results in more than 50% of patients with Helicobacter pylori infection and severe atrophic gastritis of the stomach mucosa (Kokkola et al., 1998, 2002). If the patient is on PPI medication, the 13C urea breath test gives a false negative result in almost 40% of cases. Bismuth-containing medicines and antibiotics further increase the number of false results. In these cases, too, a serological antibody test gives more reliable results than the 13C urea breath test (Graham et al. 2002, 2003, Gatta et al. 2004). The 13C urea breath test alone does not detect atrophic gastritis, and in addition, the false negative result means that the Helicobacter pylori infection and therefore also the atrophic gastritis are left untreated. These problems do not occur when the Helicobacter pylori test is performed using a serological antibody test (Kokkola et al. 1998, 2002), particularly when GastroPanel, in addition to its serological IgG&IgA H. pylori antibody test, includes biomarkers indicating the presence of atrophic gastritis in the stomach mucosa (pepsinogen I and the ratio of pepsinogen I and II: atrophic gastritis in the corpus; gastrin-17: atrophic gastritis in the antrum).
4. Kokkola A., Rautelin H, Puolakkainen P, Sipponen P, Färkkilä M, Kosunen TU. Diagnosis of Helicobacter pylori- infection in Patients with Atrophic Gastritis: Comparison of Histology, 13C Urea Breath Test, and serology. Scnd J Gastroenterol 2000;25:138-141. - ”Conclusions: H. pylori infection would

have been missed in most patients with atrophic gastritis without the analysis of *H. pylori* antibodies. Therefore, in patients with atrophic gastritis, the use of serology is encouraged in diagnosing *H. pylori* infection.”

5. Kokkola A, Rautelin H, Puolakkainen P, Sipponen P, Färkkilä M, Haapiainen R, Kosunen TU. Positive result in serology indicates active *Helicobacter pylori* infection in patients with atrophic gastritis. *J Clin Microbiol.* 1998; 36 (6) 1808-10. -“In patients with atrophic corpus gastritis, positive serology results may indicate an ongoing infection in spite of negative 13C-UBT and histology results.”
6. Sipponen P, Ranta P, Helske T, et al. Serum Levels of Amidated Gastrin-17 and Pepsinogen I in Atrophic Gastritis: An Observation Case-Control Study, *Scand J Gastroenterol* 2002;37:785-791
7. Whiting JL, Sigurdsson A, Rowlands DC, et al. The long term results of endoscopic surveillance of premalignant lesions. *Gut* 2002;50:378-81.
8. Malfertheiner P, Megraud F, O’Morain C, et al. Current concepts in the management of *Helicobacter pylori* infection – The Maastricht 2-2000 Consensus Report. *Aliment Pharmacol Ther* 2002; 16:167-180
9. Zagari RM, Nicolini G, Casanova S, et al. Diagnosis of atrophic gastritis in the general population based upon a combination of three non invasive tests, *Gut* 2002; 51 (suppl 11): A39.
10. Väänänen H, Vauhkonen M, Helske T, et al. Non-Endoscopic Diagnosis of Atrophic Gastritis with a Blood Test. Correlation between Gastric Histology and Serum Levels of Gastrin-17 and Pepsinogen I. A Multicenter Study. *Eur J Gastroenterol Hepatol* 2003; 15: 885-891
11. Pasechnikov VD, Chukov SZ, Kotelevets SM, et al. Invasive and non-invasive diagnosis of *Helicobacter pylori*-associated atrophic gastritis: A comparative study, *Scand J Gastroenterol* 2005; 40: 297-301
12. Karvonen AL, Ylävatsavaivan diagnostiikkaa: pitääkö aina tähytää ?, *Suomen Lääkärelehti* 44/2005;60:4487 – 4490.
13. Sipponen P, Vauhkonen M, Helske T, et al. Patients with Barrett’s esophagus show low circulating levels of gastrin-17, *World J Gastroenterol* 2005; 11(38): 5988-5992
14. Sipponen P, Laxen F, Huotari K, et al. Prevalence of Low Vitamin B12 and High Homocysteine in Serum in an Elderly Male Population: Association with Atrophic Gastritis and *Helicobacter pylori* infection, *Scand J Gastroenterol* 2003; 38:1209-1216
15. Salaspuro M. Ovatko “testaa ja hoida”-strategian haitat hyötyä suuremmat, *Duodecim* 2005;8:852-853
16. Salaspuro M. Dyspepsian kaksipuolisen seulontamenetelmän markkinointi on ennenaikaista, *Duodecim* 2005;121:1191-3
17. Järvinen L. Tapausselostus: GastroPanel – uusi ase dyspepsian hoidon kehittämiseen, *Yksityislääkäri* 2005; 2: 94 – 98
18. Correa P. Human gastric carcinogenesis: a multistep and multifactorial process. First American Cancer Society Award Lecture on Cancer Epidemiology and Prevention. *Cancer Res* 1992;52:6735-6740[Abstract]. - Kts. myös [www.google.com /search: Correa’s cascade](http://www.google.com/search:Correa's%20cascade).

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Di Mario F, Franze A, Cavallaro LG. Non-Invasive Diagnosis for Gastric Diseases. *One Global Medicine s.r.l.*2004; 1-48, [www.biohit.com / Literature / Dignostics](http://www.biohit.com/Literature/Diagnostics); 2004 Books

DiMario F, Cavallaro LG, Liatopoulou A, ym. Accuracy of “serological gastric biopsy” in a cohort dyspeptic patients, Poster presentation at the DDW 2005, May 15-18, in Chigago, IL, USA

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Gatta L, Perna F, Ricci C, ym. Effect of proton pump inhibitors and antacid therapy on 13C urea breath test and stool test for *Helicobacter pylori* infection. *Am J Gastroenterol* 2004;99:823-829

Graham KS, Graham DY. Contemporary Diagnosis and Management of *H. pylori* – Associated Gastrointestinal Diseases, Published by Handbooks in Health Care Co, Newtown, Pennsylvania, USA, 2002

Graham DY, Opekun AR, Hammoud F, Yamaoka Y, Reddy R, Osato MS, El-Zimaity HM. Studies regarding the mechanism of false negative urea breath tests with proton pump inhibitors. *Am J Gastroenterol.* 2003;98:1005-9.

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Manes G, Menchise A, deNucci C. Empirical prescribing for dyspepsia: a randomised controlled trial of test and treat versus omeprazole treatment. *BMJ* 2003;326:1118 – 1123

Nurgalieva Z, El-Zimaity H, Graham D, ym. Gastric atrophy in North America: Histology vs. Non-invasive testing, Poster presentation at the DDW 2005, May 15-18, in Chigago, IL, USA

Pasechnikov VD, Chukov SZ, Kotelevets SM, ym. Invasive and non-invasive diagnosis of *Helicobacter pylori*-associated atrophic gastritis: A comparative study, *Scand J Gastroenterol* 2005; 40: 297-301

Rugge M, Correa P, Dixon MF. ym. Gastric mucosal atrophy: interobserver consistency using new criteria for classification and grading. *Aliment Pharmacol Ther* 2002;16:1-12

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Sipponen P, Laxen F, Huotari K, ym. Prevalence of Low Vitamin B12 and High Homocysteine in Serum in an Elderly Male Population: Association with Atrophic Gastritis and *Helicobacter pylori* infection, *Scand J Gastroenterol* 2003; 12: 1209 –

Sipponen P, Vauhkonen M, Helske T, ym. Patients with Barrett’s esophagus show low circulating levels of gastrin-17, *World Gastroenterol* 2005, in press

Talley NJ, Vakil N, Delaney G, ym. Management issues in dyspepsia: current consensus and controversies. *Scand J Gastroenterol* 2004; 39 (10): 913-918

Uemura N, Okamoto S, Yamamoto S, ym. *Helicobacter pylori* infection and the development of gastric cancer, *N Eng J Med* 2001; 345:784-789

Varis K, Sipponen P, Laxen F ym. the Helsinki Gastritis Study Group, Implications of serum pepsinogen I in early endoscopic diagnosis of gastric cancer and dysplasia, *Scand J Gastroenterol* 2000; 9: 950-956

Väänänen H, Vauhkonen M, Helske T, ym. Non-Endoscopic Diagnosis of Atrophic Gastritis with a Blood Test. Correlation between Gastric Histology and Serum Levels of Gastrin-17 and Pepsinogen I. A Multicenter Study. *Eur J Gastroenterol Hepatol* 2003; 15: 885-891

Waldum HL, Brenna E, Sandvik AK. Long-term safety of proton pump inhibitors: risks of gastric neoplasia and infections. *Expert Opin Drug Saf* 2002;1:29-38

Whiting JL, Sigurdsson A, Rowlands DC, ym. The long term results of endoscopic surveillance of premalignant lesions. *Gut* 2002;50:378-81.

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Zagari RM, Nicolini G, Casanova S, ym. Diagnosis of atrophic gastritis in the general population based upon a combination of three non invasive tests, *Gut* 2002; 51 (suppl 11): A39.