

**BIOHIT HealthCare**

Innovating for Health

# PEPSINOGEN I

ELISA kit for the measurement of human pepsinogen I  
in EDTA plasma and serum

## Instructions for use

**REF** 601010.01

**IVD**



For *in vitro* diagnostic use  
Store at 2-8 °C Upon Receipt

Biohit Oyj Laippatie 1, FI-00880 Helsinki, Finland

Tel. +358 9 773 861, [info@biohit.fi](mailto:info@biohit.fi), [www.biohithealthcare.com](http://www.biohithealthcare.com)

Pepsinogen I REF 601010.01

CONTENTS

1. INTENDED USE.....	4
2. CLINICAL BACKGROUND.....	4
3. PRINCIPLE OF THE TEST .....	5
4. WARNINGS AND PRECAUTIONS .....	5
5. SPECIMEN COLLECTION AND HANDLING.....	6
6. KIT CONTENTS, REAGENT PREPARATION AND STABILITY FOR MATERIALS PROVIDED.....	6
6.1. Microplate .....	6
6.2. Washing Buffer Concentrate (10 x).....	7
6.3. Diluent Buffer.....	7
6.4. Blank Solution.....	7
6.5. Calibrators .....	7
6.6. Control .....	7
6.7. Conjugate Solution .....	8
6.8. Substrate Solution .....	8
6.9. Stop Solution .....	8
6.10. Incubation Covers .....	8
6.11. Instructions for Use .....	8
7. MATERIALS REQUIRED BUT NOT PROVIDED .....	8
8. STORAGE AND STABILITY .....	9
9. TEST PROCEDURE .....	9
10. RESULTS .....	11
10.1. Quality Control Values .....	11
10.2. Calculation of the Results.....	11
10.3. Prevalence.....	12
10.4. Interpretation of the Results .....	12
11. LIMITATIONS OF THE PROCEDURE.....	13
12. PERFORMANCE CHARACTERISTICS .....	13
13. REFERENCES.....	16
14. DATE OF ISSUE .....	18
15. WARRANTY .....	18
16. ORDERING INFORMATION .....	18
17. EXPLANATION OF THE SYMBOLS USED IN LABELS.....	19
18. SHORT OUTLINE OF THE PROCEDURE .....	20

APPENDIX: QUALITY CONTROL CERTIFICATE

## 1. INTENDED USE

Pepsinogen I (PGI) test is an *in vitro* microplate-based quantitative enzyme-linked immunosorbent assay (ELISA) for the determination of human pepsinogen I from EDTA plasma or serum samples. The test aids in diagnosis of an advanced atrophic gastritis in the gastric corpus from the patients at risk to develop malignant cellular changes in stomach mucosa. In addition, it aids in screening conditions that necessitate additional examination or treatment from healthy stomach mucosa.

Pepsinogen I test can also be used together with the Pepsinogen II test (REF 601020.02) by which the PGI/PGII ratio is determined. The test can be conducted either manually or automatically and is to be used by healthcare professionals.

## 2. CLINICAL BACKGROUND

This test is intended to identify patients who have an advanced atrophic gastritis in the gastric corpus and who, correspondingly, are at increased risk for gastric cancer (1, 2). The serum or plasma PGI (P-PGI, S-PGI) assay is a reliable tool for detecting patients with advanced atrophic corpus gastritis (3-6); the sensitivity and specificity of the test are 92% and 90%, respectively.

Pepsinogen I (PGI) is a precursor enzyme of pepsin and is synthesized by the chief cells and neck cells of the gastric corpus (from so-called oxyntic glands of the gastric mucosa). The major part of PGI is secreted into the gastric lumen but a small amount can be found in the blood. The P/S-PGI level reliably correlates with the number of chief cells in the gastric corpus mucosa. Correspondingly, the loss of chief cells results in a linear decrease in P/S-PGI. The loss of chief cells is, on the other hand, a result of atrophic gastritis.

For unknown reasons, atrophic gastritis increases the risk of gastric cancer, the risk being even 5-fold in patients with advanced atrophic gastritis in the corpus and even 90-fold in advanced atrophic pangastritis (both antrum and corpus affected) compared to the cancer risk in persons with normal gastric mucosa (2).

The screening of middle-aged (50-69 years), smoking men in Finland with the S-PGI test has revealed that a low S-PGI level (<25 µg/l) is detected in 9.8% of men of whom 4.7% revealed either a gastric cancer

or precancerous lesion by endoscopy (7). Corresponding results have also been published in earlier studies (8-17).

## 3. PRINCIPLE OF THE TEST

This PGI ELISA is based on a sandwich enzyme immunoassay technique with a PGI specific capture antibody adsorbed on a microplate and a detection antibody labeled with horseradish peroxidase (HRP).

The assay proceeds according to the following reactions:

1. A monoclonal antibody, specific to human PGI, on the polystyrene surface of the wells binds PGI molecules present in the sample.
2. Wells are washed to remove the residual sample.
3. An HRP-conjugated monoclonal detection antibody is added to the wells and it binds to the PGI molecules.
4. The wells are washed and TMB-substrate is added. The substrate is oxygenized by the enzyme and a blue colored end product is produced.
5. The enzyme reaction is terminated with stop solution. The solution in the microwells should turn yellow. The intensity of the yellowish color developed is directly related to the PGI concentration of the sample.

## 4. WARNINGS AND PRECAUTIONS

### For *in vitro* diagnostic use

CAUTION: Handle plasma and serum samples as potential biohazardous material.

All human blood and control samples are to be treated as potentially infectious and handled according to standard precautions (e.g., GLP, GMMP, CLSI M29). Please refer to internationally or nationally recognized manuals concerning biosafety issues, such as Laboratory Biosafety Manual by World Health Organization or Biosafety in Microbiological and Biomedical Laboratories by Centers for Disease Control and Prevention/ National Institutes of Health.

This kit contains reagents manufactured from animal blood or tissue components. All recommended precautions for the handling of bio-derivatives should be followed.

Always use protective gloves and clothing when handling patient samples. Use a safety pipetting device for all liquid transfers. Read all instructions prior to performing this assay.

Components containing ProClin may cause an allergic skin reaction (see Safety Data Sheet). Dispose of solutions containing ProClin according to local waste management legislation.

Any serious incident that occurs in relation to the use of this kit shall be reported immediately to the manufacturer (contact details in chapter 16) and the competent authority.

## 5. SPECIMEN COLLECTION AND HANDLING

Fasting for 4 hours is recommended prior to blood sampling. Blood sample is collected by venipuncture into e.g. a plastic EDTA or serum tube without additives. Plasma blood tubes are mixed immediately by turning them upside down 5-6 times and tubes for serum allowed to clot (for minimum 30 minutes) at room temperature (20...25°C). Serum after clotting and plasma immediately is separated by centrifugation (e.g. plastic tube, acceleration up to 2000 G, 10-15 minutes). Samples can be stored for 7 days in a refrigerator at 2...8°C, and 3 days at room temperature. For longer storage the samples should be stored frozen (preferably at -70°C, alternatively at -20°C). Mix the samples thoroughly after thawing. Avoid repeated freezing and thawing of the samples. Grossly hemolysed, lipemic or turbid specimens should be avoided.

Please refer to Gastrin-17 Advanced ELISA instructions for use if testing the Biohit GastroPanel ELISA assays (Pepsinogen I, Pepsinogen II, Gastrin-17 Advanced, *Helicobacter pylori* IgG antibodies) from the same sample.

## 6. KIT CONTENTS, REAGENT PREPARATION AND STABILITY FOR MATERIALS PROVIDED

The reagents are sufficient for 96 wells and three separate runs. Reagents of different kit lots should not be mixed.

### 6.1. Microplate

**Contents:** 12 x 8 strips in frame coated with high-affinity, monoclonal anti-human -PGI IgG<sub>1</sub>.

**Preparation:** Ready for use. Do not combine strips/wells from different microplates, even if they would have the same lot number.

**Stability:** Stable until expiry date. Discard the strips after use. Note that

after opening the microplate foil, some crystals may be formed in the bottom of the wells and may affect the results. Discard the crystalized strips. Do not combine strips/wells from different microplates, even if they would have the same lot number.

### 6.2. Washing Buffer Concentrate (10 x)

**Contents:** 120 ml of 10 x phosphate buffer saline (PBS) concentrate containing Tween 20 and 0.1 % ProClin 300 as preservative.

**Preparation:** Dilute 1 to 10 (e.g. 100 ml+ 900 ml) with distilled water and mix well.

**Stability:** The diluted solution is stable for two weeks refrigerated (2...8 °C).

### 6.3. Diluent Buffer

**Contents:** 100 ml of phosphate buffer containing bovine serum albumin, Tween 20, 0.1% ProClin 300 as preservative and red dye extract.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.4. Blank Solution

**Contents:** One vial containing 1.5 ml of human serum-based phosphate buffer with 0.1% ProClin 300 as preservative.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.5. Calibrators

**Contents:** Three vials each containing 1.5 ml of human serum-based calibrators with 0.1% ProClin 300 as preservative. The calibrators have lot-specific PGI values of approximately 25, 100 and 200 µg/l. The exact PGI concentration of the calibrators is labelled on the vials.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.6. Control

**Contents:** One vial containing 1.5 ml of human serum-based PGI control with 0.1% ProClin 300 as preservative. The expected PGI level of the control is indicated on the label.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.7. Conjugate Solution

**Contents:** 15 ml of HRP-conjugated monoclonal anti-human-PGI in stabilizing buffer with 0.02% methylisothiazolone, 0.02% bromonitrodioxne and 0.002% other active isothiazolones as preservatives.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.8. Substrate Solution

**Contents:** 15 ml of tetramethylbenzidine (TMB) in aqueous solution.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date. Avoid exposure to direct light.

### 6.9. Stop Solution

**Contents:** 15 ml of 0.1 mol/l sulphuric acid.

**Preparation:** Ready for use.

**Stability:** Stable until expiry date.

### 6.10. Incubation Covers

Three plastic sheets to cover the microplate during incubation.

### 6.11. Instructions for Use

## 7. MATERIALS REQUIRED BUT NOT PROVIDED

- ▶ Distilled or deionized water
- ▶ Micropipettes and disposable tips, to accurately deliver 20 - 1000  $\mu$ l
- ▶ Pipettes to accurately deliver 1-10 ml
- ▶ 8-channel pipette delivering 100  $\mu$ l
- ▶ Graduated cylinder, 1000 ml
- ▶ Vortex mixer for sample dilutions
- ▶ Test tubes for specimen dilutions
- ▶ Microplate washer
- ▶ Paper towels or absorbent paper
- ▶ Timer
- ▶ Incubator, 37 °C

- ▶ Microplate reader, 450 nm
- ▶ E. g. plastic blood collection tube for plasma or serum
- ▶ Container for ice-water bath

## 8. STORAGE AND STABILITY

Store the Pepsinogen I kit refrigerated (2...8°C). When stored at these temperatures the kit is stable until the expiration date printed on the box label and the label of each individual kit component. Do not freeze or expose the kit to high temperatures or store at above 8°C when not in use. The substrate solution is sensitive to light. The microplate or individual strips should not be removed from the foil pouch until equilibrated to room temperature (20...25°C). Unused strips must be returned to the foil pouch, sealed and stored at 2...8°C.

Do not use reagents after the expiration date printed on the label. Do not use reagents from kits with different lot numbers or substitute reagents from kits of other manufacturers. Use only distilled or deionized water. The components of the kit are provided at precise concentrations. Further dilution or other alterations of the reagents may cause incorrect results.

### Indication of Kit Deterioration

Liquid components should not be visibly cloudy or contain precipitated material. At 2...8°C, the washing buffer concentrate may, however, partially crystallize, but the crystals will dissolve by mixing at room temperature (20...25°C). The substrate solution should be colourless or pale blue/yellow. Any other color indicates deterioration of the substrate solution.

## 9. TEST PROCEDURE

### PRELIMINARY PREPARATIONS

Allow all reagents and the microplate to reach room temperature (20...25°C). Warm the incubator to 37°C. Dilute the washing buffer concentrate 1 to 10 (e.g. 100 ml + 900 ml) with distilled or deionized water. **Read the complete assay procedure before starting. It is recommended that all calibrators and samples are applied on the plate as duplicates. It is necessary to use calibrators and the control in each test run.**

Mix all the reagents and samples well before use.

#### STEP 1: SPECIMEN DILUTION

Dilute serum or plasma samples 1 to 10 (50 µl + 450 µl) with the diluent buffer, **mix well**.

#### STEP 2: SAMPLE

Mix and pipette 100 µl of the blank solution (BS), the calibrators (CAL1-CAL3), the control and diluted samples (S1, S2 etc.) into the wells as duplicates (see Figure 1). Cover the plate with the incubation cover. Incubate for 60 minutes at 37°C.

**Figure 1. Pipetting Order.**

	1	2	3	4	5
A	BS	BS			
B	CAL1	CAL1			
C	CAL2	CAL2			
D	CAL3	CAL3			
E	Control	Control			
F	S1	S1			
G	S2	S2			
H	etc.	etc.			

#### STEP 3: WASHING

Wash the wells three times with 350 µl of the diluted (1 to 10) washing buffer and gently tap the inverted plate a few times on a clean paper towel.

#### STEP 4: CONJUGATE

Pipette 100 µl of the mixed conjugate solution into the wells, preferably with an 8-channel pipette. Cover the plate with the incubation cover. Incubate for 30 minutes at 37°C.

#### STEP 5: WASHING

Wash the wells three times with 350 µl of the diluted (1 to 10) washing buffer and gently tap the inverted plate a few times on a clean paper towel.

#### STEP 6: SUBSTRATE

Pipette 100 µl of the mixed substrate solution into wells with an 8-channel pipette. Start the incubation time after pipetting the substrate solution into the first strip and continue the incubation for 30 minutes at room temperature (20...25°C). Avoid direct exposure to light during incubation.

#### STEP 7: REACTION STOP

Pipette 100 µl of the mixed stop solution with an 8-channel pipette into the wells.

#### STEP 8: MEASURING OF RESULTS

Measure the absorbance at 450 nm within 30 minutes.

## 10. RESULTS

### 10.1. Quality Control Values

Good Laboratory Practice requires the use of appropriate controls to establish that all the reagents and protocols are performing as designated. The Pepsinogen I ELISA is provided with the control serum. Quality control charts should be maintained to follow the performance of the control or appropriate statistical methods should be used for analyzing control values, which should fall within the appropriate confidence intervals employed in each laboratory.

### 10.2. Calculation of the Results

Assay results can be analyzed by using a) manual method or b) automated methods, where the absorbance readings are converted to pepsinogen I concentrations. Since the calibrators are ready to use, the results of the patient samples are not multiplied by the dilution factor.


#### a) Manual Method

Calculate the mean absorbance of the duplicate determinations of the blank solution, the calibrators, the control and samples. Subtract the mean of the blank solution from itself (consider this as the first point of the calibration curve), the calibrators, the control and samples. Graph the calibrator curve by plotting the mean absorbance for the first point and each calibrator (y-axis) against the PGI concentrations given for the calibrators (x-axis). Draw a best fit curve to construct a calibration curve. Use the mean absorbance value for each sample and the control to interpolate the PGI value from the calibration curve.

#### b) Automated Methods

There are several computer programs available for interpolating the unknown concentrations, automatically. A simple 2nd order polynomial fit is adequate for interpolating unknown concentrations within the calibrator

range. However, if sample absorbance value exceeds the absorbance value of the highest calibrator, a more complex extrapolating algorithm may be more appropriate. A typical calibration curve is shown in Figure 2.



**Figure 2. Example of a Typical Calibration Curve.**

### 10.3. Prevalence

In an elderly population (age above 50 years), 10% shows advanced atrophic corpus gastritis and abnormal S-PGI levels (S-PGI < 25 µg/l). Approximately 5% of these patients show gastric cancer or precancerous lesions in endoscopy (7).

### 10.4. Interpretation of the Results

- A low P/S-PGI result (P/S-PGI < 30 µg/l) indicates advanced (moderate and severe) atrophic gastritis of the corpus mucosa. This cut-off level has been determined using the Biohit Pepsinogen I ELISA kit based on large clinical material. A low P/S-PGI is an indicator for upper gastrointestinal endoscopy (gastroscopy) because of the increased risk of these patients of developing cancer prelesions and gastric cancer.
- It is recommended that the given limits are considered as guidelines. Also the PGI results determined for given specimen with assays from different manufacturers can vary due to differences in standardization, assay methods and reagent specificity. Results obtained from other manufacturers' assay method should not be used interchangeably.

- This Pepsinogen I ELISA assay enables wide range measurements of both low and high concentrations of P/S-PGI.

## 11. LIMITATIONS OF THE PROCEDURE

As with any diagnostic procedure the Biohit Pepsinogen I ELISA test results must be interpreted together with the patient's clinical presentation and any other information available to the physician. Samples suspected of having PGI concentrations greater than the highest calibrator should be further diluted (final dilution 1 to 20) before assay.

## 12. PERFORMANCE CHARACTERISTICS

### Within-Assay Imprecision

The within-assay imprecision was determined with four serum samples. These samples were run as 17 replicates in one run.

Sample	Mean PGI (µg/l)	CV%
1	9.3	5.9
2	31.3	3.0
3	78.1	2.4
4	137.2	2.4

### Between-Assay Imprecision

The between-assay imprecision was evaluated in six assays using four serum samples. The PGI concentration of these samples was measured as duplicates.

Sample	Mean PGI (µg/l)	CV%
1	10.3	4.9
2	47.4	3.0
3	80.4	3.4
4	133.9	2.9

### Specificity/Cross-Reactivity

The cross-reaction and interference by PGII was tested by spiking five serum samples with PGII (Company Z) at the concentrations of up to 200 µg/l. The test showed no significant increase or reduction in the signal of the samples with a PGII concentration of 200 µg/l.

Sample	PGII added (µg/l in sample)	PGI observed (µg/l)	Added/0-sample (%)
1	-	4.2	-
	10	4.2	100
	50	4.4	105
	200	4.4	105
2	-	26.3	-
	10	24.6	93.5
	50	27.3	104.0
	200	24.5	93.2
3	-	56.2	-
	10	59.6	106.0
	50	60.9	108.0
	200	53.3	94.8
4	-	72.7	-
	10	70.8	97.4
	50	72.5	99.7
	200	71.9	98.9
5	-	100.1	-
	10	98.3	98.2
	50	99.9	99.8
	200	98.7	98.6

### Sensitivity

The sensitivity of the test was determined in two different ways:

- 1) Dilution of a kit calibrator with PGI concentration 10.0 µg/l with diluent buffer gives a 10% CV limit at a concentration of 0.6 µg/l.
- 2) Mean of 25 zero replicates + 2 standard deviations corresponds to 1.9 µg/l PGI.

### Recovery

Four serum samples were spiked with 6.2, 31.0 and 60.8 µg/l human pepsinogen I (purified human PGI, Biohit Diagnostics).


The average recovery was:

- 6.2 µg/l 97.5%
- 31.0 µg/l 89.1%
- 60.8 µg/l 72.6%

### Correlation

Correlation was shown with the relationship between the serum levels of pepsinogen I and histological status of the corpus mucosa (18).

### Serum Pepsinogen I and Corpus Mucosa - Finnish Case-Control study



### Linearity

Three serum samples were assayed in serial dilutions with the diluent buffer to determine the linearity of Biohit Pepsinogen I ELISA. Results are listed in the following table.

Sample	Dilution factor	Observed (µg/l)	Expected (µg/l)	Recovery (%)
1	1	212.4	-	-
	2	105.4	106.2	99
	4	52.8	53.1	99
2	1	146.4	-	-
	2	72.8	73.2	99
	4	39.3	36.6	107
3	1	59.0	-	-
	2	29.7	29.5	101
	4	14.8	14.7	101



### 13. REFERENCES

1. Varis K. Surveillance of pernicious anemia. In *Precancerous Lesions of the Gastrointestinal Tract*. Scherlock P, Morson PC, Barbara L, Veronesi U (eds), Raven Press, New York 1983:189-194.
2. Sipponen P, Kekki M, Haapakoski J, Ihamäki T, Siurala M. Gastric cancer risk in chronic atrophic gastritis: statistical calculations of cross sectional data. *Int J Cancer* 1985; 35:173-177.
3. Varis K, Samloff IM, Ihamäki T, Siurala M. An appraisal of tests for severe atrophic gastritis in relatives of patients with pernicious anemia. *Dig Dis Sci* 1979; 24:187-191.
4. Varis K, Kekki M, Härkönen M, Sipponen P, Samloff IM. Serum pepsinogen I and serum gastrin in screening of atrophic pangastritis with high risk of gastric cancer. *Scand J Gastroenterol* 1991; 186:117-123.
5. Knight T, Wyatt J, Wilson A, Greaves S, Newell D, Hengels K, Corlett M, Webb B, Forman D, Elder J. *Helicobacter pylori* gastritis and serum pepsinogen levels in a healthy population: development of a biomarker strategy for gastric atrophy in high groups. *Br J Cancer* 1996; 73:819-824.
6. Borch K, Axelsson CK, Halgreen H, Damkjaer Nielsen MD, Ledin T, Szesci PB. The ratio of pepsinogen A to pepsinogen C: a sensitive test for atrophic gastritis. *Scand J Gastroenterol* 1989; 24:870-876.
7. Varis K, Sipponen P, Laxén F, Samloff IM, Huttunen JK, Taylor PR, Heinonen OP, Albanes D, Sande N, Virtamo V, Härkönen M, Helsinki Gastritis Study Group. Implications of serum pepsinogen I in early endoscopic diagnosis of gastric cancer and dysplasia. *Scand J Gastroenterol* 2000; 35:950-956.
8. Miki K, Ichinose M, Ishikawa KB, Yahagi N, Matsushima M, Kakei N, Tsukada S, Kido M, Ishihama S, Shimizu Y, Suzuki T, Kurokawa K. Clinical application of serum pepsinogen I and II levels for mass screening to detect gastric cancer. *Jpn J Cancer Res* 1993; 84:1086-1090.
9. Hattori Y, Tashiro H, Kawamoto T, Kodama Y. Sensitivity and specificity of mass screening for gastric cancer using the measurement of serum pepsinogens. *Jpn J Cancer Res* 1995; 86:1210-1215.
10. Yoshihara M, Sumii K, Haruma K, Kiyohira K, Hattori N, Tanka S, Kajiyama G, Shigenobu T. The usefulness of gastric mass screening using serum pepsinogen levels compared with photofluorography. *Hiroshima J Med Sci* 1997; 46:81-86.
11. Miki K, Ichinose M, Shimizu A, Huang SC, Oka H, Furihata C, Matsushima T, Takahashi K. Serum pepsinogens as a screening test of extensive chronic gastritis. *Gastroenterol Jpn* 1987; 22: 33-141.
12. Kodoi A, Yoshihara M, Sumii K, Haruma K, Kajiyama G. Serum pepsinogen in screening for gastric cancer. *J Gastroenterol* 1995; 30: 452-460.
13. Aoki K, Misumi J, Kimura T, Zhao W, Xie T. Evaluation of cut-off levels for screening of gastric cancer using serum pepsinogens and distribution of levels of serum pepsinogen I, II and of PG1/PGII ratios in a gastric cancer case-control study. *J Epidemiol* 1997; 7:143-151.
14. Kikuchi S, Wada O, Miki K, Nakajima T, Nishi T, Kobayashi O, Inaba Y. Serum pepsinogen as a new marker for gastric carcinoma among young adults. Research group on prevention of gastric carcinoma among young adults. *Cancer* 1994; 73:2695-2702.
15. Yoshihara M, Sumii K, Haruma K, Kiyohira K, Hattori N, Kitada Y, Komoto K, Tanka S, Kajiyama G. Correlation of ratio of serum pepsinogen I and II with prevalence of gastric cancer and adenoma in Japanese subjects. *Am J Gastroenterol* 1998; 93:1090-1096.
16. Farinati F, Di Mario F, Plebani M, Cielo R, Fanton MC, Valiante F, Masiero M, DeBoni M, Della Libera G, Burlin A. Pepsinogen A/ Pepsinogen C or Pepsinogen A multiplied by gastrin in the diagnosis of gastric cancer. *Ital J Gastroenterol* 1991; 23:194-206.
17. Nomura AM, Stemmermann GN, Samloff IM. Serum pepsinogen I as a predictor of stomach cancer. *Ann Intern Med* 1980; 93:537-540.
18. Sipponen P, Ranta P, Helske T, Kääriäinen I, Mäki T, Linnala A, Suovaniemi O, Alanko A, Härkönen M. Serum levels of amidated gastrin-17 and pepsinogen I in atrophic gastritis: an observational case-control study. *Scand J Gastroenterol* 2002; 37:785-791.

#### 14. DATE OF ISSUE

Pepsinogen I kit insert.  
Version 09, 05-2022

#### 15. WARRANTY

The Manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer. In case of interpretation disputes the English text applies.

This Biohit diagnostic kit has been manufactured according to our ISO 9001 / ISO 13485 quality management protocols and has passed all relevant Quality Assurance procedures related to this product.

In case of any serious incident in relation to the product, contact the manufacturer.

#### 16. ORDERING INFORMATION

Pepsinogen I ELISA test kit.  
Cat. No. 601010.01.



**BIOHIT OYJ**

Laippatie 1  
00880 Helsinki, Finland  
Tel: +358-9-773 861  
E-mail: [info@biohit.fi](mailto:info@biohit.fi)  
[www.biohithealthcare.com](http://www.biohithealthcare.com)

#### 17. EXPLANATION OF THE SYMBOLS USED IN LABELS

	CE-Marking
	For <i>in vitro</i> diagnostic use
	Catalogue number
	Batch code
	Use by
	Consult instructions for use
	Temperature limitation. Store at +2...8°C
	96 determinations
	Do not re-use
	Manufacturer

## 18. SHORT OUTLINE OF THE PROCEDURE

Allow all the reagents to reach room temperature (20...25°C)

**Remember to mix all the reagents and samples well just before pipetting**

\*

After mixing, pipette 100 µl of the blank solution, the calibrators, the control and diluted (1 to 10) patient samples into the wells

\*

Incubate for **60 min at 37°C**

\*

Wash the wells 3 times with 350 µl of the diluted washing buffer

\*

Pipette 100 µl of the mixed conjugate solution into the wells

\*

Incubate for **30 min at 37°C**

\*

Wash the wells 3 times with 350 µl of the diluted washing buffer

\*

Pipette 100 µl of the mixed substrate solution into the wells

\*

Incubate for **30 min** at room temperature (**20...25°C**)

\*

Pipette 100 µl of the mixed stop solution into the wells

\*

Read at **450 nm** within 30 minutes