# BIOHIT HealthCare

## Innovating for Health

### ENGLISH

**BIOHIT Celiac quick test for the detection of human anti-tissue transglutaminase antibodies in blood**
REF 602070 (20 tests), REF 602080 (50 tests)

**1. INTENDED USE**
BIOHIT Celiac quick test is a qualitative *in vitro* immunochromatographic test aiding in diagnosis of celiac disease by detecting antibodies (IgA/IgG/IgM) against human tissue transglutaminase in whole blood, plasma or serum samples. The test is conducted manually or semi-automatically and is intended to be used by healthcare professionals, either in a laboratory-based or point of care (POC) setting.

**2. SUMMARY AND EXPLANATION**
Celiac Disease (CD) is an enteropathy caused by a permanent intolerance to gluten and specifically to gliadin, one of its protein fragments. The ingestion of this protein in people with the genetic predisposition induces a severe injury of the intestinal mucosa that is histologically characterized by hyperplasia of the crypts, with total or subtotal atrophy of the intestinal microvilli. Though the definitive diagnosis of celiac disease is based on characteristic histological changes observed *in vivo*, serological tests such as the detection of anti-gliadin, anti-tTG and anti-endomysium antibodies, offer cost-effective and less invasive methods of disease detection. Transglutaminase is the principal auto-antigen recognised by the antiendomysial antibodies. It is particularly useful for pediatric patients (up to 16 years-old).

**3. PRINCIPLE OF THE QUICK TEST**
Anti-tTG antibodies present in a blood sample react with latex particles with human recombinant tTG. These latex particle-tTG-Anti-tTG complexes react the reaction zone through a chromatographic process, where immobilised human tissue transglutaminase captures the complex, forming a red/pink line.

**4. WARNINGS AND PRECAUTIONS**
For *in vitro* diagnostic use. Always use protective gloves when handling patient samples. Read all instructions prior to performing the quick test. Do not use test plates after their expiry date. Discard the used test plates to biohazardous waste according to local and national regulations. **CAUTION:** Handle whole blood specimens as potentially biohazardous material. All blood specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the U.S. Department of Health and Human Services (Bethesda, MD, USA) publication *Biosafety in Microbiological and Biomedical Laboratories*, 2009, 5th ed. (CDC/NIH) and No. (CDC) 21-1112 on reports of laboratory safety procedures on different diseases or any other local or national regulations.

**5. KIT CONTENTS AND REAGENT PREPARATION**
Each Celiac Quick Test kit contains reagents for 20/50 tests. Store the kit components in a refrigerator (2–8°C) until use.
5.1. Test Devices: 20/50 individually packed test devices, ready for use. Stable until the expiration date. Discard the used cassettes after use.
5.2. Diluent buffer: Diluent buffer, 5ml (REF 602070), 2\*5 ml (REF 602080), ready for use. Stable until the expiration date.
5.3. Lancets: 20/50 lancets, ready for use. Stable until the expiration date. Discard the lancets after use.
5.4. Micropipettes: 20/50 micropipettes (20 ul), ready for use. Stable until the expiration date. Discard the micropipettes after use.
5.5. Instructions for Use

**6. MATERIALS REQUIRED, BUT NOT PROVIDED**
Timer, Alcohol wipes, Gloves

**7. STORAGE AND STABILITY**
Store the Celiac Quick Test kit in a refrigerator (2–8°C). When stored at these temperatures the kit is stable until the expiration date printed on the box label and the label of test plate bags. Do not freeze or expose the kit to high temperatures or store at temperatures above 8°C when not in use. Do not use test cassettes, diluent buffer, lancets or micropipettes after the expiration date printed on the label.

**8. TEST PROCEDURE**
Read the complete assay procedure before starting.
**PRELIMINARY PREPARATION**
Bring a test cassette and the diluent buffer to room temperature. Remove the test cassette from its sealed pouch immediately before use and place on a flat surface.

**SAMPLES**
The test has been designed for use with fresh, non-hemolyzed samples of whole blood. Samples can be kept in the refrigerator for one day before testing. The test can also be used with plasma or serum samples. Prepare a 1/10 dilution of the serum or plasma in the dilution buffer. For example: 15 ul of serum/plasma + 135 ul of buffer. Mix well by gentle pipetting. Carefully add 125 ul of the dilution, drop by drop, into the center of the circular window.

**STEP 1:** Sample collection and application
**ATTENTION:** Collection and application of the sample must be swift, in order to avoid sample coagulation in the micropipette.

A whole blood sample is collected by finger prick, the index, middle or ring finger can be used. Clean the fingertip with alcohol and allow to air dry. Position the hand palm-side up and with a NEW, sterile lancet, prick the finger tip. To operate the lancet (Figure 1.1, a) Turn at least 1.5 turns of the grey lancet cap until you feel a separation from the device. This activates the lancet. Do not pull the cap – just rotate it. b) Press the open end of the lancet on the patient's fingertip and press the release button. Discard the lancet immediately into an appropriate container.

Hold the finger lower than the elbow and gently apply intermittent pressure to the base of the pricked finger several times. Wipe away the first drop of blood (using a sterile gauze pad or cotton).

To collect the blood sample (Figure 2.)
(1) Hold the blood collection micropipette horizontally, with the air vent free, (2) Touch the blood sample with the tip of the tube. Do not touch or scrape the skin. Fill the micropipette to the fill line, avoiding air bubbles. Do not squeeze the micropipette during collection. Test the sample immediately to avoid coagulation of the blood.
(3) Slowly dispense the blood sample onto the test cassette. Carefully, place the pipette into the round window (indicated by an arrow), then cover the pipette ventilation hole with a finger and gently press the upper part of the pipette and squeeze the bulb of the micropipette.
**STEP 2:** Dilution buffer (4) Wait 30–60 seconds until the blood has been absorbed and add two drops of dilution buffer in the same window. Add one drop at a time, holding the dropper bottle in a vertical position. Avoid contamination of the dropper bottle with the sample. If the dropper is not used, add 100 microliters of dilution buffer.
**STEP 3:** Reading of the results (Figure 3.) Visually read the results 10 minutes after adding the dilution buffer. Using a Casette Reader Instrument the results must be read exactly at 10 minutes for both positive and negative results (refer to the Quick Test Reader (REF 740400) Instructions).

**9. RESULTS**
The test results are to be interpreted as follows (Figure 3.):
NEGATIVE: Only one BLUE band appears across the result window close to the “C” letter (Control Line) of the test cassette. This band must always appear.
POSITIVE: In addition to the BLUE control band, a distinguishable PINK-RED band also appears across the result window close to the letter “T” (Test line) of the test cassette. The intensity of the line depends on the concentration of antibodies in the sample.

If no BLUE band appears in the control area, the test is INVALID and should be repeated with a new test cassette. Any line or colour that appears after 10 minutes has no diagnostic value. The final diagnosis should not be based only on the result obtained with one test; it should be based on a correlation of test results with other appropriate data and with clinical symptoms

**10. LIMITATIONS**

The Celiac Quick Test is qualitative, and no quantitative interpretation on the magnitude of the result should be made with respect to the intensity of a positive line. As with any diagnostic procedure, the Celiac Quick Test result must be interpreted together with the patient's clinical presentation and any other information available to the physician. It is very important to add the correct quantity of sample. Insufficient sample volume may not reach the reaction area, and the test may not run correctly. If the sample volume is too large, the reaction mix in the device will be diluted, which may result in a false negative result. It is very important to control the reaction time. If the reaction time is shorter than recommended, this may lead to false negative result. If the reaction time is longer than recommended, this may lead to a false result.

**11. WARRANTY**
Biohit shall remedy all defects discovered in any Product (the “Defective Product”) that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Biohit's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.
The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by Biohit.This Biohit Diagnostic kit has been manufactured according to ISO 9001 / ISO 13485 quality management protocols and has passed all relevant Quality Assurance procedures related to this product. In case of interpretation disputes the English text applies.
In case of any serious incident in relation to the product, contact the manufacturer.

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