BIOHIT ColonView

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2. INTENDED USE AND BACKGROUND

BIOHIT ColonView quick test is a qualitative *in vitro* immunochromatographic test aiding in diagnosis of lower gastrointestinal disorders, such as colorectal cancers and large adenomas that bleed, by detecting human hemoglobin (Hb) and hemoglobin/haptoglobin complex (Hb/Hp) in stool samples. The test is conducted manually or semi-automated and is intended to be used by healthcare professionals only (1,2).

ColonView Hb and Hb/Hp Fecal Occult Blood Test (FIT) test was designed to aid diagnosis of lower GI disorders, such as colorectal cancers and large adenomas that bleed. In 2012, colorectal cancer was the third most common cancer in men (746,000 cases) and the second in women (614,000 cases), with more than 690 000 annual deaths worldwide (3). Screening for colorectal cancer reduces disease-specific mortality by increasing the cancer detection at its' early stages (4).

3. PRINCIPLE OF COLONVIEW FIT

The Test Cassette strip is coated with anti-human hemoglobin and anti-human haptoglobin antibodies on the Test region (T) and antibodies on the Control region (C). The test result is evaluated by the intensity of the developing test bands, which are created by concentration of colloidal gold conjugate antibody to the C and T regions of the test cassettes (4).

4. SAMPLE COLLECTION AND HANDLING

Collect a random sample of stool with the Stool Collection Paper provided with the Sample Collection Package (or with a clean dry container). For the use of the Stool Collection Paper, please refer to the Sample Collection Instruction provided in the Kit. Also using Sample collection tube, please fill it as instructed. Too much stool in the stool collection tube results in an invalid result. The specimen(s) can be stored in the refrigerator (2–8 °C) for no more than 11 days, or at room temperature (max. 25°C) for no more than 5 days.

5. TEST PROCEDURE

1. The Test Cassette and the Sample Collection Tube containing the stool sample should be brought to room temperature (20...30°C) at least 10 minutes before testing.

2. Take the required number of Test Cassettes from the foil packaging only immediately before performing the test.

Mark the Test Cassette with the name of the patient or with another form of identification.

3. Carefully shake the Sample Collection Tube to ensure that the stool sample mixes properly with the sample buffer.

4. Open the white cap of the Sample Collection Tube. Take a paper towel and break the seal of the sample collection tube. The paper towel prevents the solution from splashing while opening the tube. Hold the collection tube upright and add 3 drops of the solution into both round sample windows (S) of the Test Cassette

5. Time to results is 15 minutes. Please note when reading the test lines (T), that strongly positive results may be interpreted visually (seen by eye) sooner than 15 minutes. Visually negative results must be interpreted exactly at 15 minutes. Using a Cassette Reader Instrument the results must be read exactly at 15 minutes for both positive and negative results (refer to the Quick Test Reader (REF 740400) Instructions). Results are only valid when control lines (C) are visible.

6. STORAGE AND STABILITY

The kit can be stored at 2 ... 30°C until the expiration date printed on the outside of the box.

6. QUALITY CONTROL OR INTERNAL PROCEDURAL CONTROL

The test contains a procedural control. A colored line that appears in the control region (C) shows that each test is performed correctly. The background may become slightly yellowish in color during testing, depending on the color of the stool sample.

7. TEST RESULTS AND VALID TEST

Positive or negative results when reading with Reader Instrument: Measure the Test Cassettes with the Reader Instrument according to Reader Instrument Instructions. The result is shown in the Instrument display as "positive" or "negative" after about 15 seconds.

8. WARRANTY

The kit shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in Biohit's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL. The period of this warranty is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by Biohit. This Biohit Diagnostic kit has been manufactured according to ISO 9001 / ISO 13485 quality management protocols.

In case of interpretation disputes the English text applies.

In case of any serious incident in relation to the product, contact the manufacturer.

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