BIOHIT Helicobacter pylori UFT300

for the detection of *H. pylori* in biopsy specimens and aiding of diagnosis of *H. pylori* infection.
REF 602019 (50 tests), REF 602021 (100 tests),
REF 602005PLA (5 tests with plate), REF 602019PLA (50 tests with plate).

1. INTENDED USE AND BACKGROUND

BIOHIT *Helicobacter pylori* UFT300 quick test is a qualitative *in vitro* test aiding in diagnosis of *H. pylori* infection from human gastric mucosa by detecting bacterial urease activity in gastric biopsy specimens with biochemical reaction. The test is conducted manually and is intended to be used by healthcare professionals in Near Patient Testing (NPT) settings.

Helicobacter pylori infection is the most important cause of chronic gastritis. A gastric mucosa colonized by *H. pylori* is invariably inflamed; this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life (1). The chronic inflammatory process can lead to atrophic gastritis, which has been linked with peptic ulcer and gastric cancer, two of the most important diseases of the upper gastrointestinal tract (2-5).

2. PRINCIPLE OF THE QUICK TEST

H. pylori is characterized by its capacity to produce urease, an enzyme able to break down urea into ammonium and bicarbonate, with the reaction: $(\text{NH}_4)_2\text{CO} + 2\text{H}_2\text{O} + \text{H}^+ \rightarrow 2\text{NH}_3 + \text{HCO}_3^-$. To selectively identify *H. pylori* with high sensitivity, the test exploits the pH variation caused by the production of ammonium, revealed by the color change of the pH indicator from yellow to magenta.

3. WARNINGS AND PRECAUTIONS

- All biopsy specimens should be regarded as potentially infectious and handled according to standard precautions (e.g. GLP, GMMP, CLSI M29).
- Use protective clothing. Contact of the test solution with skin or eyes should be prevented. If the solution contacts the skin or eye, flush with large amount of water. See the Safety Data Sheet for further information.
- Do not use beyond the expiration date printed on the outside of the box.
- Possible disinfection of the bottle is done by spraying or wiping the bottle surface with disinfectant whilst the bottle cap is closed. Wipe off any residual disinfectant as it may interfere with the result if in contact with the reagent.
- Dispose materials according to local waste management legislation.

4. SPECIMEN COLLECTION AND HANDLING

Patients should not have taken antibiotics or bismuth salts within three weeks, or H2-blockers and PPI drugs within two weeks, preceding endoscopy. Incomplete eradication of *H. pylori* can give false negative results. It is recommended that biopsy specimens are tested immediately after collection. Several biopsies from one patient collected from antrum and corpus can be analysed simultaneously in a test when appropriate. Do not let the specimen dry before testing. If there is visible blood or bile on the specimens, dab the biopsy briefly on sterile gauze. Even a small amount of blood can mask the color of the reaction.

Do not contaminate the biopsy forceps or biopsy specimen with lidocaine hydrochloride or formalin, or other lubricants or tissue fixatives. It is recommended that the biopsies for BIOHIT *Helicobacter pylori* UFT300 test are collected first, before any biopsies are collected for histology. If the forces become contaminated with formalin or other tissue fixatives rinse thoroughly in water before collecting the biopsies for BIOHIT *Helicobacter pylori* UFT300 test. Contamination with lidocaine hydrochloride (e.g. Xylocaine®) can inhibit the reaction and cause false negative results.

5. KIT CONTENTS AND REAGENT PREPARATION

- Test Tubes / Test Plates
Test tubes: REF 602019 50 pcs, REF 602021 100 pcs.
Test plates: REF 602005PLA 5 pcs, REF 602019PLA 50 pcs.
- Reagent Solution Bottle
Ready for use. Stable until expire date printed on the label or for 6 months when in use.
Reagent volume: REF 602005PLA 2.5 ml; REF 602019 & REF 602019PLA 6.5 ml; REF 602021 13 ml.
- Instructions for Use

6. MATERIALS REQUIRED BUT NOT PROVIDED

Forceps, timer, sterile gauze pads, gloves. Performance of the kit may be tested with BIOHIT *Helicobacter pylori* Control (REF 602017).

7. STORAGE AND STABILITY

Store the kit at room temperature (18...25 °C). The in-use stability is 6 months. Do not use after the expiration date printed on the label.

8. TEST PROCEDURE

- Write the ID of the patient on the test tube's label (provided on separate sheets)/ on plate label
- Open the tube/label covering the well of the plate and add one or more biopsies from the same patient into the tube/well.
- Add one drop of the reagent solution into the test tube/well - just the minimum amount required to cover the biopsy (Figure 1). Keep the bottle upright while dropping (Figure 2). Only if the biopsy is not fully covered by the solution, add one more drop; avoid adding too much of the reagent solution, not to dilute the reaction. Put the cap back on the tube/put the label back to cover the well of the plate.
- Shake the tube/plate by swirling it horizontally on a flat surface for 5 seconds. Make sure the biopsy and the solution are at the bottom of the test tube/well or of the plate by giving the tube/plate a light tap.

9. RESULTS

The color reaction is easy to interpret by removing the cap/label of the plate and looking into the tube from the top. The presence of *H. pylori* in the biopsy makes the color of the reagent solution turn from yellow to magenta. The presence of blood or alkaline substances may induce a variation in the yellow color without causing a turn into magenta. The color should be

read at 5 min after the addition of the reagent on the top of biopsy for both positive and negative results. Do not read the color after 5 min. Discard the tube/plate after 5 minutes reading. **Negative result:** Color of test solution remains yellow. Figure 3.

Positive result: Color of test solution turns from yellow to magenta. Figure 4.
Invalid result: If the color of the test solution turns to another color than magenta immediately after the addition of the reagent, the result can be influenced by the presence of interfering substances in the biopsy.

The time needed for the color change is inversely proportional to the *H. pylori* concentration

10. LIMITATIONS

The result of the test - as for all diagnostic tests - must be considered together with other clinical data available to the physician. False negative results can occur when the concentration of *H. pylori* in the biopsy is below the sensitivity limit of the test, or when the biopsy does not contain *H. pylori*, because of the uneven distribution of bacteria in the gastric mucosa.

Helicobacter pylori infection is the most important cause of chronic gastritis. A gastric mucosa colonized by *H. pylori* is invariably inflamed; this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life (1). The chronic inflammatory process can lead to atrophic gastritis, which has been linked with peptic ulcer and gastric cancer, two of the most important diseases of the upper gastrointestinal tract (2-5).

11. WARRANTY

Biohit shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in Biohit's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL. The period of this warranty is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by Biohit. This Biohit Diagnostic kit has been manufactured according to ISO 9001 / ISO 13485 quality management protocols. In case of interpretation disputes the English text applies. In case of any serious incident in relation to the product, contact the manufacturer.

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