


# BIOHIT *Helicobacter pylori* UFT300 quick test



ULTRA-FAST, READY-TO-USE TEST PLATE AND TUBE VERSION



## Ultra-fast *H. pylori* detection from a biopsy

The BIOHIT *H. pylori* UFT300 is a true quick test for the detection of *H. pylori* from a biopsy. The biopsy taken during gastroscopy can be tested immediately to diagnose *H. pylori* infection or to check the success of eradication therapy. The result is ready in only 5 minutes enabling diagnosis and reporting in endoscopy. The UFT300 quick test has excellent sensitivity and specificity which makes it a highly reliable tool for diagnostics.

- Ready-to-use test kit
- Results ready in 5 min
- Storage at room temperature
- Sensitivity 94,5 %, Specificity 100 %<sup>1</sup>
- Test and report immediately

## Product specifications

Sample: gastric biopsy specimen

Test principle: *H. pylori* urease activity

Kit components:

- test tube/plate and reagent
- instruction for use

Storage: 18-25 °C

Quality control: positive control available



Not available in the United States. In Japan for research use only.  
Check availability for other markets.

## Biopsy testing could not be easier

The BIOHIT *H. pylori* UFT300 quick test is easy and effortless to use. The biopsy specimen is placed into the test tube/well of the plate and mixed with the test reagent. Color change indicates the presence of *H. pylori* in the specimen. The procedure is safe and hygienic for the user, and the interpretation of the result doesn't require specialist training.

## Ordering details:



### BIOHIT *Helicobacter pylori* UFT300 quick test with plate

REF	Item
602005PLA	BIOHIT <i>Helicobacter pylori</i> UFT300, 5 tests
602019PLA	BIOHIT <i>Helicobacter pylori</i> UFT300, 50 tests


### BIOHIT *Helicobacter pylori* UFT300 quick test with tube

REF	Item
602019	BIOHIT <i>Helicobacter pylori</i> UFT300, 50 tests
602021	BIOHIT <i>Helicobacter pylori</i> UFT300, 100 tests


REF	Item
602017	Positive control reagent

## Simple test procedure with plate

1. Place 1 or more biopsies into the well of the plate.



2. Add one drop of the reagent solution. Cover the plate with label and shake for 5 seconds.




3. Read after 5 minutes.

Negative




Positive




## Simple test procedure with tube

1. Place 1 or more biopsies into the tube.



2. Add one drop of the reagent solution. Cap the tube and shake for 5 seconds.



3. Read after 5 minutes.



Negative

Positive



The result is ready in only 5 minutes.

### CONTACT

**Biohit Oyj** Laippatie 1,  
00880 Helsinki, Finland  
Tel. +358 9 773 861  
info@biohit.fi

**Biohit HealthCare Srl.** Via Figino 20/22,  
20016 Pero (Milano), Italy  
Tel +39 02 381951, Fax +39 02 38236521,  
info.italy@biohit.fi

**Biohit HealthCare Ltd.** Pioneer House,  
North Rd, Ellesmere Port, CH65 1AD,  
United Kingdom  
Tel. +44 151 550 4 550, Fax. +44 151 550 4 551  
info@biohithealthcare.co.uk