# BIOHIT

Innovating for Health

## GastroPanel<sup>®</sup> quick test NT v2.0

INSTRUCTIONS FOR USE  
INSTRUCCIONES DE USO

GastroPanel<sup>®</sup>

602410 (30 tests / pruebas)

602420 (30 tests / pruebas)

For *in vitro* diagnostic use  
Store at 2 to 30°C upon receipt

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## EXPLANATION OF THE SYMBOLS USED IN LABELS

| Symbol                                                                              | English                                |
|-------------------------------------------------------------------------------------|----------------------------------------|
|    | For <i>in vitro</i> diagnostic use     |
|    | Catalogue number                       |
|    | Batch code                             |
|    | Use by (yyyy-mm-dd)                    |
|    | Consult instructions for use           |
|    | Storage limitation. Store at 2 to 30°C |
|    | 30 determinations                      |
|    | Do not reuse                           |
|    | Keep dry                               |
|    | Keep away from heat                    |
|   | CE Mark                                |
|  | Exclamation mark, GHS07                |
|  | Manufacturer                           |
|  | Non-corrugated fiberboard (paperboard) |
|  | Unique Device Identifier               |

## INSTRUCTIONS FOR USE

English

**Note!** Other languages available at [www.biohithealthcare.com](http://www.biohithealthcare.com)

**GastroPanel® quick test NT**

**REF 602410, 602420**

## CONTENTS

|                                                                      |    |
|----------------------------------------------------------------------|----|
| 1. INTENDED USE .....                                                | 4  |
| 2. INTRODUCTION .....                                                | 4  |
| 2.1. Clinical background .....                                       | 4  |
| 2.2. Biomarkers .....                                                | 5  |
| 2.3. GastroPanel® .....                                              | 7  |
| 3. PRINCIPLE OF THE TEST .....                                       | 7  |
| 4. TRACEABILITY OF VALUES .....                                      | 8  |
| 5. WARNINGS AND PRECAUTIONS .....                                    | 8  |
| 6. KIT CONTENTS, STORAGE AND DISPOSAL OF MATERIALS PROVIDED ..       | 9  |
| 6.1. REF 602410 for EDTA-plasma and venous whole blood .....         | 9  |
| 6.2. REF 602420 for fingerprick whole blood .....                    | 10 |
| 7. MATERIALS REQUIRED BUT NOT PROVIDED .....                         | 11 |
| 8. STORAGE AND STABILITY .....                                       | 11 |
| 9. SAMPLE COLLECTION AND PREPARATION .....                           | 11 |
| 9.1. Preanalytical factors .....                                     | 11 |
| 9.2. Venous blood sampling and sample preparation .....              | 12 |
| 9.3. Fingerprick blood sampling and sample preparation .....         | 12 |
| 9.4. Gastrin-17 stimulation .....                                    | 13 |
| 10. TEST PROCEDURE .....                                             | 13 |
| 11. INTERPRETATION OF THE RESULTS .....                              | 15 |
| 11.1. General information .....                                      | 15 |
| 11.2. Reference intervals for individual biomarkers .....            | 15 |
| 11.3. Decision algorithm for GastroPanel quick test categories ..... | 16 |
| 11.4. Quality controls .....                                         | 18 |
| 12. ANALYTICAL PERFORMANCE CHARACTERISTICS .....                     | 18 |
| 13. DIAGNOSTIC PERFORMANCE .....                                     | 23 |
| 14. DATE OF ISSUE .....                                              | 24 |
| 15. WARRANTY .....                                                   | 26 |
| 16. ORDERING INFORMATION .....                                       | 26 |
| 17. REFERENCES .....                                                 | 52 |

## 1. INTENDED USE

GastroPanel quick test NT (REFs 602410, 602420) is an *in vitro* semi-automated immunological lateral flow test for the quantitative detection of pepsinogen I (PGI), pepsinogen II (PGII), gastrin-17 (G-17), and qualitative detection of antibodies against *Helicobacter pylori* (Hp), from human EDTA-plasma, and venous whole blood (REF 602410), or finger prick whole blood samples (REF 602420). GastroPanel quick test NT is used with the GP Reader NT device (REF 740450), and the embedded interpretation software. The test is to be used by healthcare professionals, either in a laboratory-based or point of care (POC) setting.

GastroPanel quick test NT is intended for diagnosing *H. pylori* infection and atrophic gastritis (AG) from patients with dyspeptic symptoms or at risk to develop malignant cellular changes in stomach mucosa. In addition, the test can aid in screening conditions that necessitate additional examination or treatment from healthy stomach mucosa.

## 2. INTRODUCTION

### 2.1. Clinical background

Dyspepsia, defined as predominant epigastric pain that has lasted at least 1 month (Moayyedi et al., 2017), has an estimated global prevalence of ca. 20% (Ford et al., 2015; Talley et al., 1992). The most commonly reported symptoms include epigastric burning, postprandial fullness, early satiety and nausea. Dyspepsia lowers the quality of life and causes economic burden both from societal and health service perspective (Maconi et al., 2002; Moayyedi & Mason, 2002). Approximately 25% of dyspepsia patients suffer from gastritis, gastroesophageal reflux disease, peptic ulcer disease or malignancy; the most important organic etiologies of dyspepsia (Harer & L. Hasler, 2019). Management of dyspepsia varies considerably in the primary health care and often a treatment is determined without a diagnostic test, or gastrointestinal endoscopy and histopathology of biopsies is used to exclude organic origin of condition (Black et al., 2018; de Jong et al., 2019). There are several implications that follow, from overuse of expensive and burdensome endoscopy to the health issues caused by a long-term use of proton pump inhibitors, or failure to diagnose and treat a *Helicobacter pylori* infection (de Jong et al., 2019; Singh et al., 2018; Sipponen & Härkönen, 2010; Zendehdel & Roham, 2019). *H. pylori* is known to be the root cause for many gastric pathologies, including gastric cancer, gastritis and ulcer, and a chronic infection caused by this pathogen may be associated with various neurological and metabolic disorders, as well as some cardiac and respiratory diseases. Currently, *H. pylori* infection is estimated to affect 50% of population globally (Holleczek et al., 2019; Zendehdel & Roham, 2019). The prevalence of gastric atrophy

is approximated to be 5-11%, varying geographically and in between age groups, increasing with age and being more common in Asian countries (Weck & Brenner, 2006). According to the Correa cascade model (Correa, 1992), gastritis is a known risk factor in the establishment of premalignant lesions that may progress to gastric cancer (de Vries et al., 2008; Koulis et al., 2019; Sipponen et al., 1985). Atrophic gastritis also predisposes to malabsorption of vitamin B12, iron, magnesium, zinc, calcium and some medicines.

## 2.2. Biomarkers

According to different anatomic regions of the stomach, gastric glands differ in morphology as well as in types of specialized cell populations. Corpus and fundus contain simple tubular glands, of which parietal cells and chief cells are responsible for gastric acid and pepsinogen I and II secretion. Antral and pyloric regions contain branched glands, which are composed of endocrine cells secreting gastrins, and mucous cells secreting pepsinogen II. Pepsinogens and gastrins are mainly excreted into the stomach lumen, but a small proportion of them diffuses into the blood stream and, hence, can be measured from a blood sample. The concentrations can be used to evaluate the condition of different regions of stomach mucosa. (Agréus et al., 2012; Miki, 2006; Syrjänen, 2016; Väänänen et al., 2003; Zagari et al., 2017).

### PEPSINOGENS I AND II

Human gastric mucosa contains aspartic proteinases that can be separated based on their physical properties into two major groups: pepsinogen I (PGI) and pepsinogen II (PGII). Pepsinogens are secreted in the stomach lumen where hydrochloric acid, secreted by the parietal cells, converts them into the corresponding active enzyme pepsins. Pepsinogen synthesis and secretion are regulated by positive and negative feedback mechanisms.

PGI is a precursor enzyme (zymogen) of pepsin, synthesized in the gastric corpus. The circulating PGI concentration closely correlates with the quantity of the chief cells in the corpus mucosa, and any loss of these cells results in decrease in the levels of PGI. As a biomarker, PGI is intended to identify patients who have mucosal atrophy (atrophic gastritis) in the gastric corpus.

Another zymogen, PGII is produced by the chief cells and mucous neck cells of the gastric corpus, in pyloric glands of the gastric antrum, and in Brunner's glands of the proximal duodenum. Hence, the level of PGII reflects the status of entire gastric mucosa. An elevated PGII level indicates mucosal inflammation, often detected in *H. pylori*-associated non-atrophic gastritis. Thus, since *H. pylori* antibody levels can remain elevated for several months even after successful eradication, PGII may be a useful marker for the confirmation of positive eradication results. PGII test complements PGI as an additional

diagnostic tool for atrophic corpus gastritis, and PGI/PGII ratio decreases with an increasing grade of atrophy.

## **GASTRIN-17**

Gastrins are linear peptide hormones produced by the G cells in the duodenum, in the pyloric part of the gastric antrum, as well as in the pancreas. The G cells release a mixture of different molecular weight gastrins into the circulation, including gastrin-71, -52, -34, -17, -14, and -6, all of which are carboxy-amidated and circulate in an O-sulfated and non-sulfated form. In healthy humans, the dominant forms of gastrin in plasma are amidated gastrin-34 (G-34) and -17 (G-17), of which G-17 is the predominant and most potent form in healthy antral tissue, and it is almost exclusively produced by the antral G cells. The main function of gastrins is to stimulate the secretion of gastric acid (HCl) by the parietal cells in the gastric corpus, as well as to increase the motility of antrum. In addition, gastrins are known to stimulate gastric chief cells to secrete pepsinogens (PGI, PGII), and induce the contraction of the lower esophageal sphincter.

The amidated form of G-17 included in this test is a direct biomarker of antral structure and function, and through a negative feedback loop, an indirect biomarker of gastric corpus. G-17 plasma levels within the normal range implicate a normal structure and function of the antrum, whereas low or high values of G-17 reflect abnormal functions of the corpus. The maximum information is obtained when G-17 testing is done separately both for fasting (G-17b) and stimulated (G-17s) levels. The measurement of plasma G-17b may also be used for the monitoring of patients who have undergone gastric surgery; secretion of G-17b is practically zero after successful radical antral resection (antrectomy). In *H. pylori* negative subjects, a low fasting level of G-17 may indicate high acid output.

## **HELICOBACTER PYLORI**

*H. pylori* is a spiral-shaped, gram-negative bacterium that colonizes in the human stomach. *H. pylori* is the most common cause of chronic infection worldwide, with a prevalence of approximately 50%; however, the majority of the infected individuals are asymptomatic. The organism is found within the mucous layer overlying the gastric epithelium, and also, within the mucosal glands, but it does not appear to invade the epithelial cells. However, the mucosa underneath and surrounding the areas of the *H. pylori* colonization is invariably inflamed; this condition is referred to as chronic superficial or non-atrophic gastritis which, if untreated, persists for life and may sequel to peptic ulceration and gastric carcinoma (de Brito et al., 2019). Antibodies are directed towards *H. pylori* and an infection can be detected and quantified from a blood sample.

## 2.3. GastroPanel®

GastroPanel was developed to meet the need to have a minimally invasive tool for identifying organic origin of dyspepsia symptoms, and to diagnose *H. pylori* infection. The levels of PGI and PGII, G-17, and antibodies to *H. pylori* to provide information on both the structure and the function of stomach mucosa, hence assisting health care professionals to treat dyspepsia patients and to screen subjects at risk to develop malignant cellular changes. (Storskrubb et al., 2008; Agréus et al., 2012; Syrjänen, 2016; Zagari et al., 2017). GastroPanel quick test NT system utilizes the same combination of validated biomarkers and decision algorithm in a Point-of-Care setting. These versions of the GastroPanel can save time and costs and expedite referral to further examinations and treatment of the patient.


Normal plasma levels of all four biomarkers indicate that the stomach mucosa has normal structure and function, whereas abnormal levels are signs of a non-healthy stomach, reflecting disturbances in the feedback mechanisms between the acid output of the corpus, PGs and G-17. For G-17 assessment, there are two options; G-17 basal (G-17b) values and G-17 stimulated (G-17s) levels, the latter being particularly important in distinguishing between functional disturbance of the antrum (high acid output) and atrophic gastritis in the antrum.

## 3. PRINCIPLE OF THE TEST

The GastroPanel quick test NT system is an immunochromatographic (lateral flow) quick test that measures the concentration of four biological markers of gastric mucosal structure and function: pepsinogen I (PGI), pepsinogen II (PGII), gastrin-17 (G-17), and antibodies to *Helicobacter pylori* (*Hp*) from EDTA-plasma, venous or fingerprick whole blood. Detection is based on time-resolved fluorescence of europium-labelled microparticles using the dedicated GP Reader NT device.

**Sample Window:**  
Sample Pad  
Conjugate Pad

**Test Window:**  
Antibody coating

The test cassette contains four parallel test strips, one for each analyte. On each test strip membrane, there are two lines with a coating. The test line contains specific antibodies to PGI, PGII, G-17, and *Hp* antigen, respectively, and the control line contains goat anti-mouse antibodies. Upon pipetting a diluted sample into the sample window, the liquid passages the sample pad and solubilizes fluorescent microparticles that are coated with specific antibodies to

PGI, PGII, G-17, or *Hp* antigen, respectively. In solution, the analytes bind to the microparticles and start migrating along the membrane by capillary action. The microparticle-bound analytes are captured by the PGI, PGII and G17 antibodies or *Hp* antigen in the test line, whereas the excess of microparticles not containing analyte is captured in the control line. The concentration of the analytes is quantified from the fluorescence intensity ratio between the test line, and the control line, and calculated by the software according to the pre-established, lot-dependent calibration by the GP Reader NT.

For installation, handling, and maintenance of the GP Reader NT, please refer to GP Reader NT User's Guide.

#### **4. TRACEABILITY OF VALUES**

The metrological traceability chain of GastroPanel quick test NT calibrators has been established according to the standard ISO 17511:2020. The standard describes methods to establish metrological traceability either to international primary measurement procedures or primary reference materials, and when these are not available, to other calibration materials and methods. As there is no primary material or method for any of the four GastroPanel quick test NT biomarkers, the reference measurement procedure and calibrators established by the manufacturer were selected as upper end of metrological traceability chain.

Values of PGI, PGII and G-17 calibrators were assigned based on the value transfer chain to previously established and validated versions of GastroPanel product family, and validated in a Clinical Comparison Study. The *H. pylori* test utilizes an in-house developed *H. pylori* antigen mixture. In GastroPanel Quick Test NT, novel calibrator and control materials were established. Their values were assigned to Biohit internal master calibrators and validated in the Clinical Comparison Study. The resulting arbitrary unit is called *Helicobacter pylori* unit, HPU.

#### **5. WARNINGS AND PRECAUTIONS**

For *in vitro* diagnostic use.

**CAUTION:** Handle blood samples as potentially biohazardous material.

All human blood and control samples are to be treated as potentially infectious and handled according to standard precautions (e.g., GLP, GMMP, CLSI M29). Please refer to internationally or nationally recognized manuals concerning biosafety issues, such as Laboratory Biosafety Manual by World Health Organization or Biosafety in Microbiological and Biomedical Laboratories by Centers for Disease Control and Prevention/National Institutes of Health.

This kit contains proteins from bovine and murine origin. Some of these might cause allergic skin reaction and/or allergy or asthma symptoms or breathing difficulties if inhaled

Always use protective gloves and clothing when handling patient samples. Use a safety pipetting device for all liquid transfers. Read all instructions prior to performing this assay.

Components containing ProClin 300 may cause an allergic skin reaction (see Safety Data Sheet). Dispose of ProClin containing solutions according to local waste management legislation.

Any serious incident that occurs in relation to the use of this kit shall be reported immediately to the manufacturer (contact details in chapter 17).

Reagents of different kit lots must not be mixed.

## 6. KIT CONTENTS, STORAGE AND DISPOSAL OF MATERIALS PROVIDED

### 6.1. REF 602410 for EDTA-plasma and venous whole blood

The kit contains reagents for 30 tests.

#### Sample Diluent Buffer

**Contents:** 4 vials with 4 ml (total 16 ml) of phosphate buffer containing casein, Tween-20, and 0.1% Proclin 300 as preservative.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date. Component may be used for 20 d after opening.

#### Test Cassette

**Contents:** 30 pieces.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date. Component must be used within 1 h after opening.

Test cassette is single use only. Dispose unused sample diluent buffer and test cassettes according to your local waste regulations. Treat used cassettes as potentially biohazardous and dispose as such according to your local waste regulations.

#### Other Items:

Quality Control Certificate - lot specific.

Instructions for use

## 6.2. REF 602420 for fingerprick whole blood

The kit contains reagents for 30 tests. Note: The kit contains 33 fingerprick sampling equipment and dilution vials, respectively, to allow possible repetition of sampling.

### Sample Diluent Buffer

**Contents:** 33 vials with 320 µl phosphate buffer containing casein, Tween-20, and 0.1% Proclin 300 as preservative.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date.

### Test Cassette

**Contents:** 30 pieces.

**Preparation:** Ready to use.

**Stability:** Stable until expiry date. Component must be used within 1 h after opening.

The sample diluent buffer vials and cassette are single use only. Dispose unused sample diluent buffer and test cassettes according to your local waste regulations. Treat used cassettes, and sample dilution vials as potentially biohazardous and dispose as such according to your local waste regulations.

## Fingerprick Sampling Equipment

| Item                    | Purpose              | Pieces in Kit |
|-------------------------|----------------------|---------------|
| Safety lancet (sterile) | Fingerprick sampling | 33            |
| 40-µl capillary pipette | Sampling device      | 66            |
| 80-µl transfer pipette  | Sample application   | 33            |

The equipment is single use only. Dispose unused items according to your local waste regulations. Treat used items as potentially biohazardous and dispose as such according to your local waste regulations.

### Other items

Quality Control Certificate – lot specific.

Instructions for use

## 7. MATERIALS REQUIRED BUT NOT PROVIDED

- GP Reader NT (REF 740450)
- Protective or disinfecting equipment for sampling
- Equipment for venous blood collection and handling [required for REF 602410]
- Centrifuge for venous blood tubes [required for REF 602410]
- Micropipette and disposable tips for accurate delivery of 80 and 320 µl [required for REF 602410]
- Empty sample dilution tube able to carry volumes up to 1.5 ml (for example, the Eppendorf tube) [required for REF 602410]
- Biohit protein powder (REF 601037 and 601038) for gastrin-17 stimulation

## 8. STORAGE AND STABILITY

The stability of the kit is established according to the results from accelerated stability studies. The stability shall be updated and notified as soon as real time stability has been established.

Store the GastroPanel quick test NT kit at 2...30°C. When stored at these temperatures, the kit is stable until the expiration date printed on the kit box label. Do not freeze or expose the kit to high humidity, or high temperatures when not in use. A test cassette should not be removed from the foil bag until equilibrated to room temperature (18...28°C). Do not use a test cassette that has been removed from the foil bag for more than 1 h ago, or has been damaged (e.g., after dropping). Do not mix test cassettes and sample diluent from kits with different lot numbers or substitute the diluent with reagents from any other source. The sample dilution tubes provided in (REF 602420) are provided in ready-to-use format. Further dilution or other alterations of the reagents may cause incorrect results.

## 9. SAMPLE COLLECTION AND PREPARATION

GastroPanel quick test NT can be used with fingerprick whole blood (REF 602420), or EDTA anticoagulated venous whole blood, and EDTA-plasma (REF 602410). Hemolyzed or grossly lipemic/cloudy samples should not be used.

### 9.1. Preanalytical factors

Some preanalytical factors may affect the measured analyte concentrations and therefore physician must be aware of these when making final diagnosis and further studies. Such are:

- Use of PPI-medication (relevance is explained in Chapter 11. Interpretation of the Results): it is recommended that PPI-medication should not be taken for 10 days prior to sampling to reach a valid interpretation of the health of stomach mucosa. However, if a pause in taking the PPI medication is not possible, physician should take this into account when making the final diagnosis.
- Fasting: it is recommended to fast at least four hours, preferably overnight, prior to sampling since the levels of gastrin-17 are very sensitive to protein ingestion. See also Chapter 9.4 for protein stimulation.

Equilibrate sample dilution buffer and test cassette to room temperature (18... 28°C) before starting the sampling and analyzing procedures.

### 9.2. Venous blood sampling and sample preparation (REF 602410)

1. Collect venous blood into an EDTA-tube (not provided) according to best practices in phlebotomy and national regulations.

Analysis of the whole blood sample must be completed within 2 hours to avoid degradation of the analytes.

**For plasma samples:** centrifuge the EDTA-tube according to manufacturer's instructions to separate plasma.

**Note:** Plasma samples can be stored frozen for analysis later. Freeze the sample immediately after separation of plasma. The sample may be stored at -20°C for two weeks, and at -70°C, for one year. Mix the samples thoroughly after thawing. The overall processing time of the sample must not exceed two hours.

2. Prefill an empty sample dilution tube (not provided) with 320 µl Diluent Buffer and transfer 80 µl of the whole blood/plasma sample into the tube (Dilution 1:5).
3. Close the tube and mix by inverting the tube gently five times.
4. Use the diluted sample immediately (Chapter 10). Do not store.

### 9.3. Fingerprick blood sampling and sample preparation (REF 602420)

1. Clean the fingertip with alcohol and allow to air dry.
2. Remove the transparent protective cap from a lancet. Position the hand palm-side up and press the open end of the lancet to the side of the fingertip to activate.
3. Wipe away the first drop of blood using a sterile gauze pad or a cotton swab.

4. Hold the finger lower than the elbow and draw blood into the 40- $\mu$ l capillary pipette (provided in the kit) keeping the pipette horizontally and gently touching the blood droplet.

Note: Hold on the pipette just under the bulb to keep the air holes on the bulb free.

Fill the sampling device up to the thin black line. Blood flow may be improved by gently applying intermittent pressure to the finger.

5. Dispense the blood **immediately** into a sample dilution tube (provided in the kit) by gently squeezing the bulb of the sampling device.
6. Repeat Steps 4. and 5. with a further 40- $\mu$ l capillary pipette to achieve the total of 80- $\mu$ l blood sample volume (Dilution 1:5). If the second capillary is not sufficiently filled, make a new prick (repeating Steps 1 to 5).
7. Close the tube and mix by inverting the tube gently five times.
8. Use the diluted sample immediately (Chapter 10). Do not store.

#### 9.4. Gastrin-17 stimulation

When a postprandial gastrin-17 analysis is needed, a protein stimulation is performed. Basal level of gastrin-17 is measured after fasting for a minimum of 4 hours (preferably overnight), after which a drink made from Biohit protein powder (REF 601037 and 601038, not provided) is consumed (see the product's own instructions for use). Postprandial level is measured 20 minutes after the ingestion. Sampling procedures are described in the previous chapters.

### 10. TEST PROCEDURE

Before using the GP Reader NT, please follow the installation instructions provided with the instrument.

1. Power up the GP Reader NT from the back of the instrument.
2. Log in and choose TEST.
3. Mark the sample number in the bottom part of the test cassette.
4. On the TEST screen, fill in the sample number (Sample No.).

By default, the calibration QR code given on the test cassette is used, i.e., the "Reagent Lot." field needs no attention.

**Note:** If the QR code on the cassette is unreadable, the batch specific QR code (found attached inside of the kit box lid) may also be used. Please refer to the user Manual of GP reader NT.

**This option should be used with caution to avoid the use of wrong batch calibration. Please ensure the correct Batch Number on the box.**

For options in employing calibration, please refer to the User Guide of GP Reader NT.

- 5.0. Depending on the applied sample type, select either “Plasma” or “Whole blood”.

For whole blood measurement, GP Reader NT applies a hematocrit correction. You may use three options:

- 5.1. If the hematocrit value of the sample is known, this may be filled in the corresponding field.

**OR**

- 5.2. Select whether a female or a male sample is to be measured. Accordingly, the GP Reader NT applies average hematocrit values of 40.2% and 45.5% for female and male samples, respectively.

**OR**

- 5.3. If the gender or hematocrit value is not known, fill in the overall average value of hematocrit, 42.85%.

**Note:** The last selection regarding plasma/whole blood measurement is retained on the instrument. Please ensure correct sample type upon subsequent measurement.

6. Select whether the analysis is performed on a fasting (basal), or a postprandial (stimulated), sample. The fasting status is important as it is used in the interpretation of test results (please refer to Chapter 11.).
7. Select PPI if the patient is using PPI-medication. This information is not used in the interpretation but is only to inform the physician.
8. Transfer 80 µl of diluted sample from sample dilution tube in each of the sample windows in the upper part of the test cassette.  
REF 602410: use a micropipette (not provided in the kit), REF 602420: use the 80-µl transfer pipette (provided in the kit).  
The thin black line indicates the 80-µl mark.
9. Incubate the test cassette at room temperature (18...28°C) for 15 minutes.

## OR

The cassette may also be incubated in the GP Reader NT: Ensure the right orientation (the arrow on the cassette pointing left) (Figure) and place the cassette on the tray.

Select “On timer” at the lower part of the screen. To activate the timer (15 min), press TEST and the reading of the cassette is automatically commenced after the indicated incubation time and you may skip Step 10. below.

10. Insert the test cassette into the GP Reader NT. Ensure the right orientation: (the arrow pointing left) (Figure) and press TEST.
11. After measurement, the results with interpretation will be displayed on the screen. The results can be printed by pressing PRINT. You can also export the data from the GP Reader NT or upload them to LIS/LIMS. See GP Reader NT instructions for more information.
12. Remove the cassette from the GP Reader NT.

## 11. INTERPRETATION OF THE RESULTS

### 11.1. General information

The interpretation should be based on all the GastroPanel markers measured from the same patient sample, and assay data must be gathered and analyzed together.

As with any diagnostic procedure, GastroPanel quick test NT results must be interpreted together with each patient's clinical presentation and any other anamnestic information available to the physician, such as as basal or stimulated sample, the use of PPI medication, and information about *H. pylori* eradication.

The GP Reader NT provides an interpretation of the results. This Test Report is shown on the screen, and it is printed together with the measured values. **Please note** that the information on proton-pump inhibitor (PPI) use is not considered in the interpretation – it is only to inform the physician. In contrast, the information on fasting status of the sample is represented in the interpretation. Chapter 11.3. provides a detailed description of the decision tree used in constituting the Test Report.

### 11.2. Reference intervals for individual biomarkers

The following reference intervals for PGI, PGII, PGI/PGII and G17 are based

on Biohit Oyj's internal data derived from GastroPanel ELISA-measurements from 7000 Finnish subjects (unpublished Biohit data). These cut-off values and reference ranges for healthy stomach were validated for GastroPanel quick test NT analytes during Clinical Comparison Study (Chapter 13).

- Pepsinogen I (PGI) 30...160 µg/l

- PGI concentration less than 30 µg/l may indicate atrophic gastritis of the corpus.

- Pepsinogen II (PGII) 3...15 µg/l

- PGII concentrations above 15 µg/l may indicate inflammation of the stomach mucosal membrane.

- Gastrin-17 fasting (G-17b) 1.8...7 pmol/l

- G-17 basal concentration after fasting less than 1.8 pmol/l may indicate high acid output.
- Concentrations after fasting above 7 pmol/l may indicate low acid output due to PPI-medication or atrophic gastritis of corpus.

- Gastrin-17 postprandial (G-17s) 3...30 pmol/l

- G-17 concentration should increase to at least 3 pmol/l after effective protein stimulation (postprandial). If it does not, it may be a sign of atrophic gastritis of the antrum.

- PGI/PGII ratio 3...20

- Ratio less than 3 may indicate atrophic gastritis of corpus.

Cut off for *Helicobacter pylori* was validated during GastroPanel quick test NT Clinical Comparison Study with 133 clinical samples.


- *Helicobacter pylori* (Hp) Antibodies < 30 HPU

- Elevated *H. pylori* antibody level may indicate an ongoing or recent infection; as antibody levels can remain elevated for several months even after successful eradication.

### 11.3. Decision algorithm for GastroPanel® quick test NT categories

In brief, the levels of the four biomarkers and one ratio are compared to each decision point, marked with black dashed rectangles dash line. The structural and functional causes of dyspeptic symptoms detected by the test system are marked with orange boxes that indicate atrophic gastritis and superficial gastritis caused by *H. pylori* infection. Green boxes represent the acid production of the stomach.

Decision tree (note: preanalytical data are not considered) PGI = pepsinogen I ( $\mu\text{g/l}$ ), PGII = pepsinogen II ( $\mu\text{g/l}$ ), PR = ratio of PGI and PGII, Hp = *Helicobacter pylori* antibodies (HPU), G17b = gastrin-17 fasting (basal) level (pmol/l), G17s = gastrin-17 stimulated (postprandial) level (pmol/l), N = no gastritis, AGA = atrophic gastritis in antrum, AGC = atrophic gastritis in corpus, AGP = atrophic gastritis in both antrum and corpus (panatropy), AGA\* = suspected atrophic gastritis in antrum, Hp+ = positive *Helicobacter pylori* finding, Hp- = *Helicobacter pylori* antibodies not detected, High Acid = Increased level of stomach acid (hyperchlorhydria), Normal Acid = Stomach acids are at the reference level, Low Acid = Decreased level of stomach acid, or no acid at all (hyperchlorhydria or achlorhydria). Low Acid may be caused by the use of proton pump inhibitors.

## 11.4. Quality controls

Good Laboratory Practice requires the use of appropriate controls to establish that all the reagents and protocols are performing as designated. Currently, the manufacturer does not have any such control available but the user is strongly recommended to establish means to monitor and control the performance. The Reagent QC feature (accessible from the Quality Control Menu in main window of the Reader) may be used to assist in implementation of controls; please refer to the Reader Manual.

If there is any concern on the performance of the test system, please contact the manufacturer.

## 12. ANALYTICAL PERFORMANCE CHARACTERISTICS

| Performance test                                                                                                   | Result                                                                                                                    |                                                                                                                             |                                                                                                                                |                  |
|--------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                                                    | Pepsinogen I                                                                                                              | Pepsinogen II                                                                                                               | Gastrin-17                                                                                                                     | <i>H. pylori</i> |
| Equivalence of specimen types                                                                                      | EDTA-plasma vs. finger prick whole blood: 20 samples at range 70 - 195 ng/ml was estimated to be 10.5% (95%CI 6.12-2.98%) | EDTA-plasma vs. finger prick whole blood: 20 samples at range 7.6 – 26.6 ng/ml was estimated to be 1.3% (95%CI -4.12-6.71%) | EDTA-plasma vs. finger prick whole blood: 20 samples at range 1.4 – 34.7 pmol/l was estimated to be -0.39% (95%CI -8.55-7.76%) | N/A              |
| Correlation with EDTA-plasma vs. finger prick whole blood and venous whole blood PGI, PGII, G-17: Bias-% (95-% CI) | EDTA-plasma vs. venous whole blood: 20 samples at range 70 - 195 ng/ml was estimated to be 11.2% (95%CI 7.8-14.6%)        | EDTA-plasma vs. venous whole blood: 20 samples at range 7.6 – 26.6 ng/ml was estimated to be -4.7% (95%CI -9.84-0.41%)      | EDTA-plasma vs. venous whole blood: 14 samples at range 1.4 – 34.7 pmol/l was estimated to be 0.7% (95%CI -5.15-6.64%)         |                  |

| Performance test                                                                                                                                                                                                                                                                                        | Result                                                                                                                                                                  |                                                                                                                                                                             |                                                                                                                                                                                     |                                                                                                                                                                                              |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                         | Pepsinogen I                                                                                                                                                            | Pepsinogen II                                                                                                                                                               | Gastrin-17                                                                                                                                                                          | <i>H. pylori</i>                                                                                                                                                                             |
| <p>Trueness</p> <p>Correlation with GastroPanel Unified ELISA (ref. 606 400) PGI, PGII, G-17: Bias-% (95-% CI), (Proportional Bland-Altman)</p> <p><i>H. pylori</i> Ab: Overall agreement (OOA95) Positive percentage agreement (PPA) (95% CI)</p> <p>Negative percentage agreement (NPA) (95-% CI)</p> | <p>160 samples at range 7.4 - 195 ng/ml was estimated to be -4.8% (95%CI -6.62- -2.98%)</p>                                                                             | <p>160 samples at range 3.5 - 54 ng/ml was estimated to be 3.6% (95%CI 1.46- 5.72%)</p>                                                                                     | <p>G17 with 160 samples at range 1.1 - 34.7 pmol/l was estimated to be 15.3% (95%CI 11.4-19.3%)</p>                                                                                 | <p>56 positives and 102 negatives, NPA was estimated to be 98% (93.1-99.5%), and PPA 92.9% (83-97.2%), and overall agreement 96.2% (93.2-99.2%)</p>                                          |
| <p>Precision</p> <p>Repeatability, within-run CV % (95% CI)</p> <p>2-way nested ANOVA, n=80</p>                                                                                                                                                                                                         | <p>20 µg/l: 8.5 (7.0-10.9) %</p> <p>36.2 µg/l: 5.6 (4.6-7.1) %</p> <p>46.1 µg/l: 5.6 (4.6-7.1) %</p> <p>94.5 µg/l: 4.1 (3.4-5.2) %</p> <p>168 µg/l: 5.7 (4.7-7.3) %</p> | <p>2.2 µg/l: 9.6 (7.9-12.3) %</p> <p>3.3 µg/l: 7.7 (6.3-9.8) %</p> <p>7.8 µg/l: 9.7 (8.0-12.4) %</p> <p>14.7 µg/l: 10.3 (8.5-13.2) %</p> <p>40.1 µg/l: 8.3 (6.8-10.6) %</p> | <p>1.9 pmol/l: 9.3 (7.7-11.9) %</p> <p>2.9 pmol/l: 6.9 (5.6-8.8) %</p> <p>6.2 pmol/l: 8.2 (6.7-10.5) %</p> <p>11.1 pmol/l: 6.9 (5.7-8.9) %</p> <p>25.2 pmol/l: 8.8 (7.2-11.3) %</p> | <p>9.89 HPU: 13.7 (11.3 – 17.5) %</p> <p>21.5 HPU: 7.89 (6.48 – 10.1) %</p> <p>45.4 HPU: 10.2 (8.41 – 13.1) %</p> <p>89.4 HPU: 15.8 (13 – 20.3) %</p> <p>127.0 HPU: 18.1 (14.9 – 23.2) %</p> |
| <p>Precision</p> <p>Repeatability, between-run CV % (95% CI)</p> <p>2-way nested ANOVA, n=80</p>                                                                                                                                                                                                        | <p>20 µg/l: 2.8 (N/A) %</p> <p>36.2 µg/l: 2.5 (1.2-40.7) %</p> <p>46.1 µg/l: 1.2 (N/A) %</p> <p>94.5 µg/l: 2.5 (1.4-9.1) %</p> <p>168 µg/l: 0 (N/A) %</p>               | <p>2.2 µg/l: 0 (N/A) %</p> <p>3.3 µg/l: 2.7 (1.1- &gt;100) %</p> <p>7.8 µg/l: 0 (N/A) %</p> <p>14.7 µg/l: 0 (N/A) %</p> <p>40.1 µg/l: 2.9 (1.2- &gt;100) %</p>              | <p>1.9 pmol/l: 0 (N/A) %</p> <p>2.9 pmol/l: 4.4 (2.6-13.9) %</p> <p>6.2 pmol/l: 2.0 (N/A) %</p> <p>11.1 pmol/l: 6.2 (4.1-12.6) %</p> <p>25.2 pmol/l: 0 (N/A) %</p>                  | <p>9.89 HPU: 10.3 (6.4 – 25.1) %</p> <p>21.5 HPU: 5.48 (3.32 – 15) %</p> <p>45.4 HPU: 8.01 (5.09 – 18.5) %</p> <p>89.4 HPU: 9.8 (5.64 – 33.8) %</p> <p>127.0 HPU: 11.2 (6.47 – 38.6) %</p>   |

| Performance test                                                                  | Result                                                                                                                                           |                                                                                                                                                                |                                                                                                                                                                    |                                                                                                                                                                       |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                   | Pepsinogen I                                                                                                                                     | Pepsinogen II                                                                                                                                                  | Gastrin-17                                                                                                                                                         | <i>H. pylori</i>                                                                                                                                                      |
| Precision<br>Repeatability, between-day CV % (95% CI)<br>2-way nested ANOVA, n=80 | 20 µg/l: 0 (N/A) %<br>36.2 µg/l: 0 (N/A) %<br>46.1 µg/l: 2.9 (1.7-8.9) %<br>94.5 µg/l: 2.9 (1.8-7.6) %<br>168 µg/l: 2.4 (1.4-8.8) %              | 2.2 µg/l: 11.6 (8.7-17.6) %<br>3.3 µg/l: 9.5 (6.9-15.4) %<br>7.8 µg/l: 8.0 (5.6-13.6) %<br>14.7 µg/l: 9.8 (7.1-15.6) %<br>40.1 µg/l: 10.8 (7.8-17.3) %         | 1.9 pmol/l: 8.4 (6.1-13.6) %<br>2.9 pmol/l: 5.0 (3.1-13.0) %<br>6.2 pmol/l: 8.7 (6.2-14.4) %<br>11.1 pmol/l: 5.8 (3.5-16.0) %<br>25.2 pmol/l: 4.0 (2.3-13.7) %     | 9.89 HPU: 0 (N/A) %<br>21.5 HPU: 9.02 (6.24 – 16.2) %<br>45.4 HPU: 6.18 (3.29 – 33.5) %<br>89.4 HPU: 0 (N/A) %<br>127.0 HPU: 0 (N/A) %                                |
| Precision<br>Repeatability, within-lab CV % (95% CI)<br>2-way nested ANOVA, n=80  | 20 µg/l: 9.0 (7.8-10.7) %<br>36.2 µg/l: 6.1 (5.3-7.2) %<br>46.1 µg/l: 6.4 (5.5-7.7) %<br>94.5 µg/l: 5.6 (4.7-6.8) %<br>168 µg/l: 6.2 (5.3-7.4) % | 2.2 µg/l: 15.1 (12.3-19.5) %<br>3.3 µg/l: 12.5 (10.2-16.1) %<br>7.8 µg/l: 12.6 (10.5-15.5) %<br>14.7 µg/l: 14.2 (11.8-17.8) %<br>40.1 µg/l: 13.9 (11.3-18.1) % | 1.9 pmol/l: 12.5 (10.5-15.7) %<br>2.9 pmol/l: 9.6 (8.1-11.8) %<br>6.2 pmol/l: 12.1 (9.9-15.3) %<br>11.1 pmol/l: 10.9 (9.2-13.6) %<br>25.2 pmol/l: 9.7 (8.3-11.6) % | 9.89 HPU: 17.1 (14.7 – 20.5) %<br>21.5 HPU: 13.2 (10.9 – 16.8) %<br>45.4 HPU: 14.4 (12.2 – 17.6) %<br>89.4 HPU: 18.6 (16 – 22.2) %<br>127.0 HPU: 21.3 (18.4 – 25.5) % |
| Precision<br>Reproducibility CV % (95% CI)<br>2-way nested ANOVA, n=80            | N/A                                                                                                                                              | N/A                                                                                                                                                            | N/A                                                                                                                                                                | N/A                                                                                                                                                                   |
| Measuring range<br>Linearity<br>(> LoQ) linearity within<br>+/- 20% non-linearity | 7.8...195 µg/l                                                                                                                                   | 1.5...54.0 µg/l                                                                                                                                                | 1.0...34.7 pmol/l                                                                                                                                                  | 10.3...210 HPU                                                                                                                                                        |
| Cross-Reactivity<br>Paired-difference testing and dose-response studies.          | Negligible 120 µg/l of PGI, and 100 pmol/l of G-17                                                                                               | Negligible to 700 µg/l of PGI, and 100 pmol/l G-17                                                                                                             | Negligible to 200 µg/l of PGI, and 50 µg/l of PGI                                                                                                                  | Negligible to presence of various micro-organisms *                                                                                                                   |

| Performance test                                            | Result                                                                                                                                   |                                                                                                                                          |                                                                                                                                          |                                                                                                                                                                    |
|-------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                             | Pepsinogen I                                                                                                                             | Pepsinogen II                                                                                                                            | Gastrin-17                                                                                                                               | <i>H. pylori</i>                                                                                                                                                   |
| Interference<br>Paired-difference testing **.               | Negligible to 1 g/dl hemoglobin, 37 mmol triglycerides, 0.3 mg/ml bilirubin, 1000 IU/ml rheumatoid factor, 3 mg/ml Na <sub>2</sub> -EDTA | Negligible to 1 g/dl hemoglobin, 37 mmol triglycerides, 0.3 mg/ml bilirubin, 1000 IU/ml rheumatoid factor, 3 mg/ml Na <sub>2</sub> -EDTA | Negligible to 1 g/dl hemoglobin, 37 mmol triglycerides, 0.3 mg/ml bilirubin, 1000 IU/ml rheumatoid factor, 3 mg/ml Na <sub>2</sub> -EDTA | Negligible to 0.775 g/dl hemoglobin, 16.9 mmol/l triglycerides, 475 umol/l conjugated bilirubin, 55.5 mmol/l glucose, 7.22 umol/l famotide, 24.3 umol/l omeprazole |
| Analytical sensitivity                                      |                                                                                                                                          |                                                                                                                                          |                                                                                                                                          |                                                                                                                                                                    |
| Limit of blank (LoB): according EP17                        | LoB: 3.3 µg/l                                                                                                                            | LoB: 0.9 µg/l                                                                                                                            | LoB: 0.7 pmol/l                                                                                                                          | LoB 6.6 HPU                                                                                                                                                        |
| Limit of detection (LoD) according to EP17                  | LoD: 4.7 µg/l                                                                                                                            | LoD: 1.5 µg/l                                                                                                                            | LoD: 1.0 pmol/l                                                                                                                          | LoD 10.3 HPU                                                                                                                                                       |
| (LoQ) according to EP17 (estimation from precision profile) | LoQ: 4.7 µg/l                                                                                                                            | LoQ: 1.5 µg/l                                                                                                                            | LoQ 1.0 pmol/l                                                                                                                           | LoQ 10.3 HPU                                                                                                                                                       |
| Room temperature variability                                | No difference in performance between 18 to 28°C                                                                                          | No difference in performance between 18 to 28°C                                                                                          | No difference in performance between 18 to 28°C                                                                                          | N/A                                                                                                                                                                |
| Humidity<br>Bias ≤ 20%                                      | No difference in performance between rH = (30...80%).                                                                                    | No difference in performance between rH = (30...80%).                                                                                    | No difference in performance between rH = (30...80%).                                                                                    | N/A                                                                                                                                                                |

\* *Helicobacter pylori*, *Campylobacter* sp., *Escherichia coli*, *Proteus* sp., *Candida albicans*, *Enterococcus faecalis*, *Staphylococcus aureus*, *Staphylococcus epidermidis*.

\*\* HAMA blockers (two specific and non-specific IgG) are used in the test and thus, HAMA interference has not been tested.

\*\*\* 300 µg/ml acetyl-salicylic acid, 500 µg/ml paracetamol, 300 µg/l indomethacin, 1 mg/ml naproxen, 100 µg/l diclofenac, 200 µg/l ibuprofen, 2 mg/l nimesulide, 20 000 U/ml penicillin.

### 13. DIAGNOSTIC PERFORMANCE

Clinical performance characteristics of GastroPanel quick test NT for atrophic gastritis were established with 135 clinical samples against a diagnostic upper-gastrointestinal endoscopy and biopsies according to Updated Sydney System (USS, Dixon et al., 1996). As the decision algorithm of GastroPanel products is optimized to detect moderate and severe atrophies, the USS-diagnostic category of mild atrophy was reduced to correspond to No atrophy-category in the correlation calculations. The decision algorithm and reference values of previous GastroPanel-product have been validated and published in numerous scientific articles (e.g. Germaná et al., 2005; Iijima et al., 2009; Storskrubb et al., 2008; Väänänen et al., 2003). The performance of this assay has not been established in paediatric populations.

The clinical performance for *Helicobacter pylori* was re-evaluated for the new Hp-antigen mix. The diagnostic accuracy of Hp with the histopathological examination of the corresponding biopsies (gold standard) as a comparative method was evaluated with 133 clinical samples to validate medical decision level i.e. cut-off value for the updated *H. pylori* test.

| Target name | TP | TN | FP | FN | Specificity<br>95% Confidence Interval | Sensitivity<br>95% Confidence Interval | PPV   | NPV   |
|-------------|----|----|----|----|----------------------------------------|----------------------------------------|-------|-------|
| AG          | 33 | 87 | 7  | 8  | 92.6%<br>(85.4 - 96.3) %               | 80.5%<br>(66.0 - 89.8) %               | 82.5% | 91.6% |
| N           | 66 | 53 | 5  | 11 | 91.4%<br>(81.4 - 96.3) %               | 85.7%<br>(76.2 - 91.8) %               | 93.0% | 82.8% |
| HP+         | 34 | 79 | 13 | 7  | 85.9%<br>(77.3-91.6) %                 | 82.9%<br>(68.7 – 91.5) %               | 72.3% | 91.9% |

TN = true negative, FP = false positive, FN = false negative, PPV = positive predictive value TP/(TP+FP) for balanced dataset, NPV = negative predictive value TN/(FN+TN) for balanced dataset.

## 14. DATE OF ISSUE

### **GastroPanel® quick test NT Instructions for Use**

Version 2.3, 2025-09-10

## 15. WARRANTY

The Manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer's specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS FOR USE.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer. In case of interpretation disputes the English text applies.

This Biohit diagnostic kit has been manufactured according to ISO 9001/ISO 13485 quality management protocols and has passed all relevant Quality Assurance procedures related to this product.

In case of any serious incident in relation to the product, contact the manufacturer.

## 16. ORDERING INFORMATION

### **GastroPanel® quick test NT, REF 602410**

(Plasma and venous whole blood, 30 tests).

### **GastroPanel® quick test NT, REF 602420**

(Fingerprick whole blood, 30 tests).

## 17. REFERENCES

See the references page 52.

## This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

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## 17. REFERENCES / REFERENCIAS

1. Agréus, L., Kuipers, E. J., Kupcinskis, L., Malfertheiner, P., Di Mario, F., Leja, M., Sung, J. (2012). Rationale in diagnosis and screening of atrophic gastritis with stomach-specific plasma biomarkers. *Scandinavian Journal of Gastroenterology*, 47(2), 136–147. <https://doi.org/10.3109/00365521.2011.645501>
2. Black, C. J., Houghton, L. A., & Ford, A. C. (2018). Insights into the evaluation and management of dyspepsia: recent developments and new guidelines. *Therapeutic Advances in Gastroenterology*, 11. <https://doi.org/10.1177/1756284818805597>
3. Correa, P. (1992). Human Gastric Carcinogenesis: A Multistep and Multifactorial Process—First American Cancer Society Award Lecture on Cancer Epidemiology and Prevention. *Cancer Research*, 52(24), 6735–6740. Retrieved from <http://cancerres.aacrjournals.org/content/52/24/6735.abstract>
4. de Brito, B. B., da Silva, F. A. F., Soares, A. S., Pereira, V. A., Santos, M. L. C., Sampaio, M. M., de Melo, F. F. (2019). Pathogenesis and clinical management of *Helicobacter pylori* gastric infection. *World Journal of Gastroenterology*, 25(37), 5578–5589. <https://doi.org/10.3748/wjg.v25.i37.5578>
5. de Jong, J. J., Lantinga, M. A., & Drenth, J. P. (2019). Prevention of overuse: A view on upper gastrointestinal endoscopy. *World Journal of Gastroenterology*, 25(2), 178–189. <https://doi.org/10.3748/wjg.v25.i2.178>
6. de Vries, A. C., van Grieken, N. C. T., Looman, C. W. N., Casparie, M. K., de Vries, E., Meijer, G. A., & Kuipers, E. J. (2008). Gastric cancer risk in patients with premalignant gastric lesions: a nationwide cohort study in the Netherlands. *Gastroenterology*, 134(4), 945–952. <https://doi.org/10.1053/j.gastro.2008.01.071>
7. Dixon, M., Genta, R., Yardley, J., & Correa, P. (1996). Classification and grading of gastritis. The updated Sydney System. International Workshop on the Histopathology of Gastritis, Houston 1994. *The American Journal of Surgical Pathology*, 20(10), 1161–1181.
8. Ford, A. C., Marwaha, A., Sood, R., & Moayyedi, P. (2015). Global prevalence of, and risk factors for, uninvestigated dyspepsia: a meta-analysis. *Gut*, 64(7), 1049 LP – 1057. <https://doi.org/10.1136/gutjnl-2014-307843>
9. Germaná, B., Di Mario, F., Cavallaro, L. G., Moussa, A. M., Lecis, P., Liatoupolou, S., Franzé, A. (2005). Clinical usefulness of serum pepsinogens I and II, gastrin-17 and anti-*Helicobacter pylori* antibodies in the management of dyspeptic patients in primary care. *Digestive and Liver Disease*, 37(7), 501–508. <https://doi.org/10.1016/j.dld.2005.01.016>
10. Harer, K., & L. Hasler, W. (2019). A Diagnostic Approach to Dyspepsia: A Clinical Casebook. [https://doi.org/10.1007/978-3-030-01117-8\\_6](https://doi.org/10.1007/978-3-030-01117-8_6)
11. Holleczeck, B., Schöttker, B., & Brenner, H. (2019). *Helicobacter pylori* infection, chronic atrophic gastritis and risk of stomach and esophagus cancer: Results from the prospective population-based ESTHER cohort study. *International Journal of Cancer*. <https://doi.org/10.1002/ijc.32610>
12. Iijima, K., Abe, Y., Kikuchi, R., Koike, T., Ohara, S., Sipponen, P., & Shimosegawa, T. (2009). Serum biomarker tests are useful in delineating between patients with gastric atrophy and normal, healthy stomach. *World Journal of Gastroenterology*, 15(7), 853–859. <https://doi.org/10.3748/wjg.15.853>
13. Koulis, A., Buckle, A., & Boussioutas, A. (2019). Premalignant lesions and gastric cancer: Current understanding. *World Journal of Gastrointestinal Oncology*, 11(9), 665–678. <https://doi.org/10.4251/wjgo.v11.i9.665>

14. Maconi, G., Tosetti, C., Stanghellini, V., Bianchi Porro, G., & Corinaldesi, R. (2002). Dyspeptic symptoms in primary care. An observational study in general practice. *European Journal of Gastroenterology & Hepatology*, 14(9), 985–990. <https://doi.org/10.1097/00042737-200209000-00009>
15. Miki, K. (2006). Gastric cancer screening using the serum pepsinogen test method. *Gastric Cancer*, 9(4), 245–253. <https://doi.org/10.1007/s10120-006-0397-0>
16. Moayyedi, P., & Mason, J. (2002). Clinical and economic consequences of dyspepsia in the community. *Gut*, 50(suppl 4), iv10 LP-iv12. [https://doi.org/10.1136/gut.50.suppl\\_4.iv10](https://doi.org/10.1136/gut.50.suppl_4.iv10)
17. Moayyedi, PM, Lacy, B., CN, A., Enns, R., CW, H., & Vakil, N. (2017). ACG and CAG Clinical Guideline: Management of Dyspepsia. *Am J Gastroenterol.*, 112(7), 988–1013. <https://doi.org/10.1038/ajg.2017.154>
18. Singh, A., Cresci, G. A., & Kirby, D. F. (2018). Proton Pump Inhibitors: Risks and Rewards and Emerging Consequences to the Gut Microbiome. *Nutrition in Clinical Practice*, 33(5), 614–624. <https://doi.org/10.1002/ncp.10181>
19. Sipponen, P., & Härkönen, M. (2010). Hypochlorhydric stomach: a risk condition for calcium malabsorption and osteoporosis? *Scandinavian Journal of Gastroenterology*, 45(2), 133–138. <https://doi.org/10.3109/00365520903434117>
20. Sipponen, P., Kekki, M., Haapakoski, J., Ihamäki, T., & Siurala, M. (1985). Gastric cancer risk in chronic atrophic gastritis: Statistical calculations of cross-sectional data. *International Journal of Cancer*, 35(2), 173–177. <https://doi.org/10.1002/ijc.2910350206>
21. Storskrubb, T., Aro, P., Ronkainen, J., Sipponen, P., Nyhlin, H., Talley, N. J., Agréus, L. (2008). Serum biomarkers provide an accurate method for diagnosis of atrophic gastritis in a general population: The Kalixanda study. *Scandinavian Journal of Gastroenterology*, 43(12), 1448–1455. <https://doi.org/10.1080/00365520802273025>
22. Syrjänen, K. (2016). A Panel of Serum Biomarkers (GastroPanel®) in Non-invasive Diagnosis of Atrophic Gastritis. Systematic Review and Meta-analysis. *Anticancer Research*, 36(10), 5133–5144. Retrieved from <http://ar.iarjournals.org/content/36/10/5133.abstract>
23. Talley, N. J., Zinsmeister, A. R., Schleck, C. D., & Melton, L. J. (1992). Dyspepsia and dyspepsia subgroups: A population-based study. *Gastroenterology*, 102(4), 1259–1268. [https://doi.org/10.1016/0016-5085\(92\)90764-P](https://doi.org/10.1016/0016-5085(92)90764-P)
24. Väänänen, H., Vauhkonen, M., Helske, T., Kääriäinen, I., Rasmussen, M., Tunturi-Hihnal, H., Sipponen, P. (2003). Non-endoscopic diagnosis of atrophic gastritis with a blood test. Correlation between gastric histology and serum levels of gastrin-17 and pepsinogen I: a multicentre study. *European Journal of Gastroenterology & Hepatology*, 15(8). Retrieved from [https://journals.lww.com/eurojgh/Fulltext/2003/08000/Non\\_endoscopic\\_diagnosis\\_of\\_atrophic\\_gastritis.9.aspx](https://journals.lww.com/eurojgh/Fulltext/2003/08000/Non_endoscopic_diagnosis_of_atrophic_gastritis.9.aspx)
25. Zagari, R. M., Rabitti, S., Greenwood, D. C., Eusebi, L. H., Vestito, A., & Bazzoli, F. (2017). Systematic review with meta-analysis: diagnostic performance of the combination of pepsinogen, gastrin-17 and anti-*Helicobacter pylori* antibodies serum assays for the diagnosis of atrophic gastritis. *Alimentary Pharmacology & Therapeutics*, 46(7), 657–667. <https://doi.org/10.1111/apt.14248>
26. Zendehdel, A., & Roham, M. (2019). Biological evidence of the relationship between *Helicobacter pylori* and associated extragastric diseases. *Journal of Cellular Biochemistry*, 120(8), 12128–12140. <https://doi.org/10.1002/jcb.28681>
27. Weck M, Brenner H. (2006). Prevalence of chronic atrophic gastritis in different parts of the world. *Cancer Epidemiol Biomarkers Prev*. 2006 Jun;15(6):1083-94. <https://doi.org/10.1158/1055-9965.EPI-05-0931>

## NOTAS

This image shows a single page of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page, leaving margins at the top and bottom. There is no handwriting or other markings on the paper.

[illegible]

**Biohit Oy**

Laippatie 1

FI-00880 Helsinki

Finland

Tel: +358 9 773 861

E-mail: [info@biohit.fi](mailto:info@biohit.fi)

[www.biohithealthcare.com](http://www.biohithealthcare.com)