

New, unique *H. pylori* ELISA test!



Helicobacter pylori IgG ELISA kit (Cat. No. 601 040.02)

The new Biohit *H. pylori* IgG antibody ELISA test finds *H. pylori* infection more reliably due to the new, unique *H. pylori* antigen.

Detected Ig Class	Sensitivity*	Specificity*	PPV*	NPV*	Clinical efficiency*
IgG	96.0%	97.4%	94.9%	98.0%	96.9%

* Against *H. pylori* Culture (N=523)

Advantages of the kit:

- Improved, unique antigen in both kits guarantees the quality of detection
- *H. pylori* infection is found more reliably
- Assay incubation time is 1.5 hours
- Shelf life is 12 months
- e-Lisa XL assay template is available



Headquarters Biohit Oyj

Laipatie 1 • 00880 Helsinki • Finland • Tel. +358-9-773 861 • Fax: +358-9-773 86 290 • info@biohit.com • www.biohit.com