

BIOHIT HealthCare

Innovating for Health

Peptest™ Cube Reader



The new Cube Reader for PepTest™ is now available!

- Quantitative results
- The cube-shaped device has an edge length of approx. 41 mm, making it probably the world's smallest reader that will even fit into your pocket, it weighs only 40 g.
- The Cube is inexpensive and very user-friendly since it has only one operating button.
- The clear test result is displayed after only **3 seconds**.

Device Specifications

The Cube is a mobile measuring device for colorimetric tests based on reflectance measurements to capture the optical density. The built-in techniques allow a quantitative and qualitative evaluation of the Peptest™ lateral flow assay. The test specific data is transferred wirelessly by radio-frequency identification (RFID) before each measurement.

The measurement data can be stored internally. Every saved result includes a unique measurement ID, test name, lot number, distributor/manufacturer, and date and time. The measurement data can be extracted via USB using a special USB cable and the Cube Data Reader software.

The device is powered by battery but can also be powered via USB.

The Cube can be used as a mobile handheld device or as a desktop measuring device remotely controlled via USB cable and the free-to-use Cube Data Reader software.



Cat. No.	Item
PEPCBRDR	PEP TEST CUBE READER

BIOHIT HealthCare

Innovating for Health

www.biohithealthcare.com

Distributed by:

Biohit HealthCare Srl
Via Boncompagni, 2
20139 Milano - Italy
Tel. +39 02 38238113
Fax +39 02 38236521
mail: info.italy@biohit.fi