

Biohit HealthCare

Innovating for Health

GastroPanel® quick test

Revolutionary quick test innovation of the unique GastroPanel® is now available!

GastroPanel® quick test is the further development of the unique Biohit GastroPanel examination. GastroPanel quick test allows fast diagnosis of and screening for *Helicobacter pylori*, atrophic gastritis with related risks as well as high acid output of stomach in symptomatic and asymptomatic patients. The GastroPanel quick test speeds up the referral to further examinations. GastroPanel® findings - i) *Helicobacter pylori* infection, ii) atrophic gastritis or iii) high acid output are indications for gastroscopy and biopsy examination.

GastroPanel quick test saves costs and unnecessary clinical appointments.

It significantly improves the patient safety, as it speeds up the referral to further examinations, treatment and follow-up.

- **Sample:** EDTA-plasma
- **Test principle:** immunological lateral-flow test for detecting GastroPanel biomarkers
- **Kit components:** test cassette, sample dilution tubes with buffer, instructions for use.

