

GastroPanel® quick test NT

Revolutionary quick test innovation of the unique GastroPanel® NT is now available!

GastroPanel® quick test is the further development of the unique Biohit GastroPanel® examination. GastroPanel quick test allows fast diagnosis of and screening for *Helicobacter pylori*, atrophic gastritis with related risks as well as high acid output of stomach in symptomatic and asymptomatic patients. The GastroPanel quick test speeds up the referral to further examinations. GastroPanel findings - i) *Helicobacter pylori* infection, ii) atrophic gastritis or iii) high acid output are indications for gastroscopy and biopsy examination.

GastroPanel® quick test saves costs and unnecessary clinical appointments.

It significantly improves the patient safety, as it speeds up the referral to further examinations, treatment and follow-up.

- **Sample:** EDTA-plasma, fingerprick whole blood
- **Test principle:** immunological lateral-flow test for detecting GastroPanel® biomarkers



GastroPanel® quick test NT		
REF	Product	Qty
602410	For plasma and venous whole blood	30 tests
602420	For fingerprick whole blood	30 tests

In the United States and Japan for research use only.

GastroPanel® quick test NT	
REF	Product
6740450	GP Reader NT Fluorescence Immunoanalyzer

In the United States and Japan for research use only.