

PRODUCT NAME AND CLONE:	Monoclonal antibody to VINCULIN , Clone FB11
HOST:	Mouse
IMMUNOGEN:	Monoclonal antibody to vinculin is derived from the hybridoma produced by fusion between myeloma cells and Balb/c spleen cells. Vinculin purified from human smooth muscle was used as an immunogen.
CLASS AND SUBCLASS:	IgG ₁
PRODUCT SPECIFICITY:	In immunoblotting the antibody identifies a 130 kD polypeptide.
PRODUCT CODE:	610014: 100 µg/ml vial 1ml 610041: 100 µg/ml vial 10ml
PRODUCT BUFFER:	PBS solution containing 1.0% (w/v) BSA and 0.09% (w/v) NaN ₃
STORAGE:	+2...+8°C
APPLICATIONS:	Immunohistochemical staining, immunoblotting and other immunoassays for cellular and tumor biology.
GUIDELINES FOR DILUTIONS:	Working dilution in western blotting is at least 1:1000.

For research use and further manufacturing only.

DS14 ver 7