

DATA SHEET

PRODUCT NAME AND CLONE:	The original BIOHIT Monoclonal antibody against GASTRIN-17 , Clone G55D4.1.
HOST:	Mouse
IMMUNOGEN:	Monoclonal antibody to gastrin-17 is derived from the hybridoma produced by fusion between myeloma cells and Balb/c spleen cells. Synthetic (modified) gastrin-17 coupled to thyroglobulin was used as an immunogen.
CLASS AND SUBCLASS:	IgG ₁
PRODUCT SPECIFICITY:	The antibody is specific to gastrin-17 with no cross reactivity to gastrin-34 or cholecystokinin.
PRODUCT CODE:	610084: 100 µg/ml vial 100 µg 610086: 100 µg/ml vial 1 mg
PRODUCT BUFFER:	PBS solution containing 1.0% (w/v) BSA and 0.09% (w/v) NaN ₃ . Other MAb concentrations and buffer compositions are available upon request. Matched Pairs: Gastrin-17 610061 and 610063.
STORAGE:	2 to 8°C
APPLICATIONS:	Immunohistochemical staining and enzyme-immunoassays.
REFERENCES:	

The product is for research use and further manufacturing only. Use in human clinical diagnostics is the responsibility of the user.

DS84 ver 1