



Product specifications:

- Detects both 25OH vitamin D2 and D3 for a clinically meaningful assessment of vitamin D status
- Calibrated against the ID-LC/MS-MS Reference Measurement Procedure
- Proven ELISA technology, support for both automated and manual analysis methods



Not available in the United States. In Japan for research use only.
Check availability for other markets.

BIOHIT Total 25OH Vitamin D



Ordering details:

REF	Product
602310.02	BIOHIT Total 25OH Vitamin D, 96 tests

Contact our team for more information.



Biohit Oyj
Laippatie 1, 00880 Helsinki, Finland
Tel. +358 9 773 861
info@biohit.fi

www.biohithealthcare.com



BIOHIT Total 25OH Vitamin D

For a conclusive determination of vitamin D status



490097en-05.2 ©BIOHIT OYJ 08/2025

BIOHIT
Innovating for Health

Multifunctional vitamin D

Vitamin D has multiple roles in the human body. In addition to its well established role in the regulation of calcium absorption and promoting bone growth, it is recognized for other health benefits including reducing risk of diseases such as type 1 diabetes and common cancers.

Clinical background

Two major forms of vitamin D are important in the human body, vitamin D2 (ergocalciferol) and vitamin D3 (cholecalciferol). Vitamin D3 is the form synthesized in the skin in response to sunlight UVB exposure yet can also be obtained from animal based foods. Vitamin D2, on the other hand, is the form synthesized by plants and is obtained from plant derived foodstuff. Both are metabolized in the liver to their respective 25OH vitamin D3 and D2 forms which are in turn converted to their active form in the kidneys.

"The best indicator of vitamin D status is the serum concentration of 25OH vitamin D."



The best indicator of vitamin D status is the serum concentration of 25OH vitamin D. For a correct diagnosis of vitamin D deficiency, insufficiency or intoxication, the assay must recognize both the D2 and D3 forms. As both forms have clinical relevance, the total concentration provides information that can be acted upon.


BIOHIT Total 25OH Vitamin D kit

The BIOHIT Total 25OH vitamin D kit is a quantitative immunoenzymatic assay. While detecting both 25OH vitamin D2 and D3, the kit provides clinically relevant information on the vitamin D status. Reliability of the results is ensured by validation against the ID-LC/MS-MS Reference Measurement Procedure (Ghent method) as approved by the Vitamin D Standardization Program (VDSP) with $R > 0.97$.

The assay protocol


The assay protocol includes just five steps making it easy to perform and reducing the hands-on time. Featuring pretreatment technology performed in assay's microplate wells, any open ELISA platform can be easily programmed to run the assay from start to finish.

1. Add 50 μ l of sample, calibrators and controls, and 150 μ l of incubation buffer to the wells.




Incubate 120 min at RT with shaking → Wash plate

2. Add 200 μ l conjugate




Incubate 30 min at RT with shaking → Wash plate

3. Add 100 μ l TMB



Incubate 15 min at RT with shaking

4. Add 100 μ l STOP



5. Read Absorbance at 450 nm